

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
DQB1 allele ^{1,2,3}																			
*04:01:01-04:01:04	1	2																	
*04:02:01-04:02:07, 04:02:09-04:02:10, 04:19, 04:21, 04:26, 04:28-04:29, 04:31-04:32, 04:34-04:35	1		3																
*04:02:08			3				7												
*04:03:01	1		3	4															
*04:03:02	1			4															
*04:04, 04:12 ⁴	1		3		5														
*04:05-04:06 ⁴	1	2			5														
*04:07	1	2										12							
*04:08	1	2						8											
*04:09	1		3						9										
*04:10	1		3							10									
*04:11	1		3								11								
*04:13	1		3										13						
*04:14	1	2							9										
*04:15	1	2								11									
*04:16	1	2			6														
*04:17	1	2												14					
*04:18	1		3									12							
*04:20	1				6														
*04:22, 04:24	1		w																
*04:23	1		3							11					15				
*04:25N	1		3												15				
*04:27	1		3													16			
*04:30			3																
*04:33	1																	17	
*04:36N	1		3																18
*03:06				4	5			?											
*03:07-03:08, 03:10:02-03:15, 03:17:01-03:18, 03:20, 03:23:01-03:23:02, 03:25:02-03:26, 03:37, 03:40, 03:48, 03:52-03:71, 03:74-03:78, 03:81-03:82, 03:101-03:112, 03:118N-03:131, 03:133-03:163, 03:165-03:167, 03:170-03:179, 03:183-03:185, 03:187-03:189, 03:192-03:194, 03:201-03:205, 03:207-03:214, 03:216-03:221, 03:223-03:230								?											
*03:19:01-03:19:02								8											
*03:22, 03:96										11					15				
*03:25:01				4	5														
*03:30, 03:72, 03:100, 06:02:01-06:02:25, 06:33, 06:47, 06:51:01, 06:106-06:107, 06:109, 06:111-06:117, 06:127, 06:130, 06:188, 06:200							7												
*03:132	1							?											
*03:215							7	?											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19
*05:01:19, 06:03:03										w									
*06:04:01-06:05:01, 06:07:01-06:07:02, 06:09:01-06:09:06, 06:12, 06:18:02, 06:21, 06:34, 06:36, 06:38, 06:41-06:42, 06:52, 06:118:01, 06:129, 06:186, 06:189											w				w				
*06:05:02, 06:84							7			w					w				

Negative Control

*06:06, 06:08:01-06:08:03, 06:11:01-06:11:03, 06:17-06:18:01, 06:25-06:28, 06:30-06:32:02, 06:35, 06:40, 06:45, 06:53:01-06:67, 06:85, 06:89, 06:91-06:94, 06:118:02, 06:120-06:121, 06:124, 06:128, 06:131-06:135, 06:137, 06:140-06:145, 06:149, 06:153-06:155, 06:157-06:158N, 06:160, 06:164-06:165, 06:167-06:172, 06:176-06:177, 06:180-06:181, 06:184, 06:190-06:191, 06:193N-06:196, 06:199, 06:201-06:204, 06:207											?									Negative Control
*06:10, 06:13:01-06:16, 06:19:01-06:20, 06:22:01-06:24, 06:29, 06:37, 06:46, 06:48-06:50, 06:51:02, 06:68-06:83, 06:95-06:97, 06:119, 06:122-06:123, 06:125-06:126, 06:136, 06:138-06:139, 06:146:01-06:147, 06:150-06:152, 06:156, 06:159, 06:161-06:163, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:192, 06:197-06:198, 06:206, 06:208						7					?									
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	