

DQ4 2019 7H9

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21
DQB1 allele ^{1,2,3}																					
*04:01:01:01-04:01:05, 04:42	1	2																		19	
*04:02:01:01, 04:02:01:04-04:02:07, 04:02:09-04:02:14, 04:19, 04:21, 04:26, 04:28-04:29, 04:31-04:32, 04:34-04:35, 04:37, 04:39-04:40	1		3																	19	
*04:02:08			3				7													19	
*04:03:01	1		3	4																19	
*04:03:02	1			4																19	
*04:03:03				4																19	
*04:04	1		3		5																
*04:05	1	2			5																
*04:06	1	2			5															19	
*04:07	1	2										12								19	
*04:08	1	2					8													19	
*04:09	1		3						9											19	
*04:10	1		3							10										19	
*04:11	1		3								11									19	
*04:12	1		3		5															19	
*04:13	1		3										13							19	
*04:14	1	2							9											19	
*04:15	1	2									11									19	
*04:16	1	2				6														19	
*04:17	1	2												14						19	
*04:18	1		3									12								19	
*04:20	1				6															19	
*04:22, 04:24	1		w																	w	
*04:23	1		3								11				15					19	
*04:25N	1		3												15					19	
*04:27	1		3													16				19	
*04:30			3																	19	
*04:33	1																		17		
*04:36N	1		3																18	19	
*04:38	1	2																		w	

Negative Control

DQ4 2019 7H9

*06:13:02, 06:15:01-06:15:02, 06:19:01, 06:22:01-06:24, 06:29, 06:50, 06:51:02, 06:69:01-06:76, 06:79:02-06:83, 06:95-06:97, 06:119, 06:122-06:123, 06:126, 06:136, 06:138-06:139, 06:146:01- 06:147, 06:150-06:152, 06:156, 06:159, 06:161-06:163, 06:166, 06:173-06:175, 06:178-06:179N, 06:182-06:183, 06:192, 06:197- 06:198, 06:206:01, 06:208, 06:211, 06:213, 06:215, 06:232, 06:235- 06:236, 06:249							7				?										
*06:84, 06:206:02							7				w										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21