

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
<b>DQB1 allele</b> <sup>1,2,3</sup>																						
*04:01:01:01-04:01:05, 04:38, 04:42, 04:61-04:62, 04:69, 04:71, 04:73N-04:74, 04:90	1	2																	19			
*04:02:01:01, 04:02:01:04-04:02:07, 04:02:09-04:02:15, 04:02:17-04:02:18, 04:02:20-04:02:21, 04:19, 04:21, 04:26, 04:28-04:29, 04:31-04:32, 04:34-04:35, 04:37, 04:39-04:40, 04:43-04:47, 04:50-04:60, 04:63-04:68N, 04:70, 04:72, 04:75-04:89	1		3																19			
*04:02:08			3				7												19			
*04:02:16, 04:30, 04:48			3																19			
*04:02:19, 04:02:22, 04:49	1																		19			
*04:03:01	1		3	4															19			
*04:03:02	1			4															19			
*04:03:03				4															19			
*04:04	1		3		5																	
*04:05	1	2			5																	
*04:06	1	2			5															19		
*04:07	1	2										12								19		
*04:08	1	2						8												19		
*04:09	1		3						9											19		
*04:10	1		3							10										19		
*04:11	1		3								11									19		
*04:12	1		3		5															19		
*04:13	1		3										13							19		
*04:14	1	2							9											19		
*04:15	1	2								11										19		
*04:16	1	2			6															19		
*04:17	1	2											14							19		
*04:18	1		3																	19	21	
*04:20	1				6															19		
*04:22, 04:24	1		w																	19		
*04:23	1		3							11					15					19		
*04:25N	1		3												15					19		
*04:27	1		3													16				19		
*04:33	1																17			19		
*04:36N	1		3																18	19		
*04:41N	1		3																	19	20	
*02:108, 02:144, 03:22:01:01-03:22:02, 05:226										11												
*03:06				4	5			?														
*03:07, 03:14:01-03:15, 03:17:01-03:18, 03:20, 03:23:01, 03:23:03, 03:25:02-03:26, 03:48, 03:52-03:54, 03:56-03:67, 03:69-03:70, 03:74-03:78, 03:81-03:82, 03:101-03:104, 03:106-03:112, 03:118N, 03:120-03:131, 03:133-03:135, 03:137, 03:140-03:149, 03:151-03:163, 03:165-03:167, 03:170-03:171, 03:173-03:179, 03:183, 03:185, 03:188-03:189, 03:192-03:194, 03:201-03:205, 03:207-03:210, 03:212-03:214, 03:216-03:221, 03:223-03:228, 03:230, 03:232-03:238, 03:240, 03:244, 03:255-03:262, 03:267-03:270, 03:273-03:275, 03:278, 03:351-03:353, 03:356N, 03:359, 03:387, 03:425, 03:433, 03:437-03:438, 03:441, 03:449-03:450, 03:455-03:456								?														

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
*03:19:01:01-03:19:05, 03:71, 03:150, 03:187, 03:263:01:01-03:263:01:02, 03:286, 03:335, 03:360, 03:376N, 03:381, 03:394, 03:401								8														
*03:25:01				4	5																	
*03:30, 03:72, 03:100, 06:02:01:01-06:02:32, 06:02:34-06:02:58, 06:10, 06:13:01, 06:14:01-06:14:03, 06:15:02-06:16, 06:19:02-06:20, 06:33, 06:37, 06:46-06:49, 06:51:01, 06:68, 06:75N, 06:77N-06:79:01, 06:106-06:107, 06:109, 06:111-06:117, 06:122, 06:125, 06:127, 06:130, 06:188, 06:200, 06:216N, 06:219, 06:224-06:228, 06:237, 06:240, 06:242, 06:255-06:256, 06:262, 06:264, 06:270:02-06:271, 06:273, 06:284, 06:286, 06:289-06:290, 06:293-06:298, 06:300, 06:304N, 06:306N, 06:308N, 06:311, 06:314-06:315, 06:317N, 06:322:01-06:322:03, 06:324, 06:326, 06:333, 06:335, 06:337-06:338, 06:341N, 06:344, 06:347, 06:354-06:357, 06:363-06:364, 06:366, 06:370, 06:372, 06:376, 06:380, 06:383N-06:384, 06:388, 06:390, 06:395, 06:397N, 06:401-06:402, 06:404-06:406, 06:408-06:409, 06:411-06:413, 06:416Q-06:417, 06:422N							7															
*03:96, 03:251											11			15								
*03:132	1																					
*03:215							7	?														
*03:395														14								
*05:01:19, 06:03:03										w												
*05:103, 05:231, 06:04:01:01-06:04:15, 06:04:17-06:05:02, 06:07:01-06:07:02, 06:09:01:01-06:09:12, 06:12, 06:18:02, 06:21, 06:34, 06:36, 06:38, 06:41:01:01-06:42, 06:52, 06:58, 06:118:01, 06:129, 06:142, 06:164, 06:171, 06:186, 06:189, 06:217, 06:254, 06:261, 06:265-06:267, 06:275, 06:280-06:283, 06:287-06:288, 06:291-06:292, 06:299, 06:301, 06:303N, 06:313, 06:318-06:319, 06:325, 06:332, 06:339, 06:343, 06:348-06:349, 06:351, 06:353, 06:358, 06:361, 06:369, 06:375, 06:381, 06:385, 06:398, 06:407, 06:414N, 06:420-06:421												w										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22

Negative Control

