

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
DQB1 allele ^{1,2,3}																						
*04:01:01:01-04:01:06, 04:38, 04:42, 04:61-04:62, 04:69, 04:71, 04:73N-04:74, 04:90, 04:92-04:95	1	2																	19			
*04:02:01:01, 04:02:01:04-04:02:07, 04:02:09-04:02:15, 04:02:17-04:02:18, 04:02:20-04:02:21, 04:02:23-04:02:24, 04:19, 04:21, 04:26, 04:28-04:29, 04:31-04:32, 04:34-04:35, 04:37, 04:39-04:40, 04:43-04:47, 04:50-04:60, 04:63-04:68N, 04:70, 04:72, 04:75-04:89, 04:91	1		3																19			
*04:02:08			3				7												19			
*04:02:16, 04:30, 04:48			3																19			
*04:02:19, 04:02:22, 04:49	1																		19			
*04:03:01	1		3	4															19			
*04:03:02	1			4															19			
*04:03:03				4															19			
*04:04	1		3		5																	
*04:05	1	2			5																	
*04:06	1	2			5														19			
*04:07	1	2										12							19			
*04:08	1	2						8											19			
*04:09	1		3						9										19			
*04:10	1		3							10									19			
*04:11	1		3								11								19			
*04:12	1		3		5														19			
*04:13	1		3										13						19			
*04:14	1	2							9										19			
*04:15	1	2									11								19			
*04:16	1	2				6													19			
*04:17	1	2												14					19			
*04:18	1		3																19		21	
*04:20	1					6													19			
*04:22, 04:24	1		w																19			
*04:23	1		3								11				15				19			
*04:25N	1		3												15				19			
*04:27	1		3													16			19			
*04:33	1																17		19			
*04:36N	1		3															18	19			
*04:41N	1		3																19	20		
*02:108, 02:144, 03:22:01:01-03:22:02, 05:226											11											
*03:06				4	5			?														
*03:07, 03:14:01-03:15, 03:17:01-03:18, 03:20, 03:23:01, 03:23:03, 03:25:02-03:26, 03:48, 03:52-03:54, 03:56-03:67, 03:69-03:70, 03:74-03:78, 03:81-03:82, 03:101-03:104, 03:106-03:112, 03:118N, 03:120-03:125, 03:127-03:131, 03:133-03:135, 03:137, 03:140-03:149, 03:152-03:163, 03:165-03:167, 03:170-03:171, 03:173-03:179, 03:183, 03:185, 03:188-03:189, 03:192-03:194, 03:201-03:205, 03:207-03:210, 03:212-03:214, 03:216-03:221, 03:223-03:228, 03:230, 03:232-03:238, 03:240, 03:244, 03:255-03:262, 03:267-03:270, 03:273-03:275, 03:278, 03:351-03:353, 03:356N, 03:359, 03:387, 03:425, 03:437-03:438, 03:441, 03:449-03:450, 03:455-03:456, 03:492								?														
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
*03:19:01:01-03:19:06, 03:71, 03:150, 03:187, 03:263:01:01-03:263:01:02, 03:286, 03:335, 03:360, 03:376N, 03:381, 03:394, 03:401, 03:485								8														
*03:25:01				4	5																	

Negative Control

[illegible]