

INTERPRETATION TABLE						
DQB1*05 SSP subtyping						
Amplification patterns of the DQB1*05 alleles						
	Well					
	1	2	3	4	5	6
Length of spec.	220	135	120	90	120	185
PCR product						
Length of int.	515	430	430	515	430	430
pos. control¹						
5'-primer²	26	26	30	30	30	38
	5'-g gg ^{3'}	5'-g gg ^{3'}	5'-gA C ^{3'}	5'-gA C ^{3'}	5'-gA T ^{3'}	5'-C gC ^{3'}
3'-primer³	87	57	57	47	57	86
	5'-g gT ^{3'}	5'-C AA ^{3'}	5'-gCT ^{3'}	5'-Cg A ^{3'}	5'-gCT ^{3'}	5'-A Cg ^{3'}
Well No.	1	2	3	4	5	6
DQB1 allele						
*050101-050102	1	2				
*050201	1		3			
*050202	1		3	4		
*050301-050302	1			4		
*0504	1				5	
*0505	1		3			6
DQB1 allele						
Well No.	1	2	3	4	5	6