

Well No.	1	2	3	4	5	6	7	8
DQB1 allele								
*05:01:01-05:01:03	1	2						
*05:02:01, 05:02:03	1		3					
*05:02:02	1		3	4				
*05:03:01-05:03:04	1			4				
*05:04	1				5			
*05:05	1		3			6		
*05:06	1			4			7	
*05:07	1	2					7	
*05:08	1			4		6		
*05:09	1			4				8
*05:10	1			4	5			
*05:11	1	2				6		
*03:01:04						6		
Length of spec. PCR product(s) ¹	225	135	120	95	120	135	180	190
					185	185		
Well No.	1	2	3	4	5	6	7	8
Specific Reaction								
1st match: DQB1*								
2nd match: DQB1*								
Control Reaction								
Failed Controls								
False Positive								
False Negative								