

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DQB1 allele ^{1,2,3}																								
*05:01:01-05:01:08	1	2																						
*05:01:09	1	2		4																				
*05:02:01, 05:02:03-05:02:07	1		3																					
*05:02:02	1		3	4																				
*05:03:01-05:03:10	1			4																				
*05:04	1				5																			
*05:05:01-05:05:02	1		3			6																		
*05:06:01	1			4			7																	
*05:07	1	2					7																	
*05:08	1			4					9															
*05:09	1			4				8																
*05:10	1			4	5																			
*05:11:01	1	2				6																		
*05:12	1	2								10														
*05:13	1			4							11													
*05:14	1		3									12												
*05:15	1			4									13											
*05:16	1			4										14										
*05:17	1		3												15									
*05:18	1	2							9															
*05:19	1		3																				23	
*05:20	1	2														16								
*05:21																		19						
*05:22	1	2																					23	
*05:23	1									10														
*05:24	1			4														18						
*05:25	1	2										12					17							
*05:26	1		3																		21			
*05:27, 05:29 ⁴	1	2					8																	
*05:28	1			4															19					
*05:30	1	2																19						
*05:31	1	2															17							
*05:32	1	2									11									20				
*05:33	1		3										13											
*05:34			3										13											
*05:35	1		3								11													
*05:36	1		3						9															
*05:37	1		3												15						21			
*05:38	1			4															20					
*05:39	1			4																	21			
*05:40	1			4									13	14										
*05:41N	1			4								12												
*05:42	1			4							11								20					
*05:43	1																					22		
*05:44		2											13											
*05:45	1	2									11													
*05:46	1		3														17							
*05:47	1		3													16								
*03:21											11													
*06:51:01, 06:66, 06:96, 06:99:01																						22		
*06:103												12												
*06:105																16								
*06:111																	17							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control