

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29
DQB1 allele ^{1,2,3}																													
*05:01:01:01-05:01:08, 05:01:11-05:01:20, 05:48, 05:51, 05:55, 05:61, 05:68, 05:70, 05:74, 05:76, 05:81, 05:89, 05:92-05:93, 05:95	1	2																											
*05:01:09-05:01:10	1	2		4																									
*05:02:01, 05:02:03-05:02:13, 05:53, 05:59, 05:65, 05:86, 05:94	1		3																										
*05:02:02	1		3	4																									
*05:03:01:01-05:03:05, 05:03:07-05:03:13, 05:03:15-05:03:16, 05:58, 05:64, 05:67, 05:85, 05:91	1			4																									
*05:03:06	1			4											15														
*05:03:14	1														15						22								
*05:04, 05:77	1				5																								
*05:05:01-05:05:02	1		3			6																							
*05:06:01-05:06:02, 05:50	1			4			7																						
*05:07	1	2					7																						
*05:08	1			4					9																				
*05:09	1			4				8																					
*05:10	1			4	5																								
*05:11:01-05:11:02	1	2				6																							
*05:12	1	2								10																			
*05:13	1			4							11																		
*05:14	1		3									12																	
*05:15	1			4									13																
*05:16	1			4										14															
*05:17	1		3												15														
*05:18	1	2							9																				
*05:19	1		3																								23		
*05:20	1	2														16													
*05:21																				19									
*05:22	1	2																									23		
*05:23	1									10																			
*05:24	1			4															18										
*05:25	1	2										12					17												
*05:26	1		3																			21							
*05:27, 05:29 ⁴	1	2						8																					
*05:28	1			4																	19								
*05:30	1	2																			19								
*05:31	1	2															17												
*05:32	1	2									11											20							
*05:33	1		3										13																
*05:34			3										13																
*05:35	1		3								11																		
*05:36	1		3						9																				
*05:37	1		3												15							21							
*05:38	1			4																	20								
*05:39	1			4																		21							
*05:40	1			4									13	14															
*05:41N	1			4								12					17												
*05:42	1			4							11										20								
*05:43	1																										22		
*05:44		2											13																
*05:45	1	2									11																		
*05:46	1		3														17												
*05:47	1		3												16														
*05:49	1	2											13																
*05:52	1				5		7																						
*05:54		2																											
*05:56	1			4																									28
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29
*05:57	1		3																	20									
*05:60				4															19										
*05:62	1	2																		20									
*05:63	1	2																			21								
*05:66	1			4		6																							
*05:69	1	2												14															
*05:71	1	2								10															24				
*05:72		2																	18										
*05:73		2												14															
*05:75	1	2																							25				
*05:78	1			4																						26			
*05:79	1		3																								27		
*05:80	1	2																						24					
*05:82																					21								
*05:83			3																										
*05:84	1	2										12																	
*05:87Q	1		3					8																					
*05:88	1	2																	18										
*05:90N	1		3									12								17									
*05:96	1			4							11	12																	
*03:03:05															15														
*03:10:02, 03:183, 06:01:12, 06:51:01, 06:66, 06:96, 06:99:01, 06:172																								22					
*03:21												11																	
*04:02:07																							w						
*04:28				4																									
*06:89																									24				
*06:103												12																	
*06:105																				16									
*06:111																				17									
*06:129						6																							
*06:146:01						6															20								
*06:156, 06:162, 06:169				4										14															
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29

Negative Control