

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
DQB1 allele ^{1,2,3}																																	
*05:01:01:01-05:01:08 , 05:01:11-05:01:38, 05:27 , 05:48, 05:55, 05:61, 05:68, 05:70, 05:74, 05:76, 05:81, 05:84, 05:89:01-05:89:02, 05:92-05:93, 05:95, 05:99, 05:114-05:115, 05:120, 05:122, 05:124, 05:126-05:127, 05:133, 05:137-05:139, 05:141, 05:144, 05:148, 05:150-05:152, 05:154-05:155, 05:157-05:158, 05:160, 05:162-05:164, 05:166-05:168, 05:171, 05:173, 05:176-05:177, 05:180, 05:183-05:185N, 05:187-05:188, 05:190, 05:193-05:195, 05:197, 05:215N-05:217, 05:219, 05:223, 05:225-05:226, 05:228, 05:230, 05:232, 05:234, 05:237																																	
*05:01:09-05:01:10	1	2		4																													
*05:02:01:01-05:02:01:06 , 05:02:03-05:02:12, 05:02:14-05:02:17, 05:02:19-05:02:22, 05:53, 05:59, 05:65, 05:86, 05:94, 05:100, 05:113, 05:118, 05:123, 05:125, 05:135-05:136, 05:142, 05:165, 05:170, 05:172, 05:174, 05:178-05:179, 05:181, 05:186, 05:189, 05:192, 05:196, 05:198-05:199, 05:222, 05:227, 05:229, 05:238-05:239, 05:241	1			3																											29		
*05:02:02, 05:02:18	1		3	4																											29		
*05:02:13	1		3																														
*05:03:01:01-05:03:05 , 05:03:07-05:03:13, 05:03:15-05:03:19, 05:03:21-05:03:23, 05:43:02, 05:58, 05:64, 05:67, 05:85, 05:91, 05:109, 05:121, 05:129, 05:140, 05:147, 05:149, 05:161, 05:191, 05:200, 05:202-05:214, 05:220-05:221, 05:224N, 05:233, 05:235N, 05:243	1			4																													
*05:03:06	1			4											15																		
*05:03:14 , 05:03:20	1														15						22												
*05:04 , 05:77, 05:146	1				5																										29		
*05:05:01-05:05:02	1		3			6																										29	
*05:06:01- 05:06:02 , 05:50	1			4			7																										
*05:07	1	2					7																										
*05:08	1			4					9																								
*05:09	1			4				8																									
*05:10	1			4	5																												
*05:11:01- 05:11:02	1	2			6																												
*05:12, 05:242 ⁴	1	2								10																							
*05:13, 05:201 ⁴	1			4							11																						
*05:14	1		3									12																			29		
*05:15	1			4									13																				
*05:16, 05:101⁴	1			4										14																			
*05:17	1		3												15																	29	
*05:18	1	2						9																									
*05:19	1		3																			23									29		
*05:20 , 05:159 ⁴	1	2													16																		
*05:21																				19													
*05:22	1	2																				23											
*05:23	1							10																									
*05:24	1			4															18														
*05:25	1	2								12							17																
*05:26	1		3																		21										29		
*05:28	1			4																													
*05:29	1	2						8																									
*05:30	1	2																		19													
*05:31	1	2															17																
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*05:32	1	2									11									20													
*05:33	1		3										13																29				
*05:34			3										13																29				
*05:35, 05:117 ⁴	1		3								11																		29				
*05:36	1		3						9																				29				
*05:37	1		3											15							21								29				
*05:38, 05:119	1			4																	20												
*05:39	1			4																		21											
*05:40	1			4									13	14																			
*05:41N	1			4								12					17																
*05:42	1			4								11									20												
*05:43:01, 05:131	1																					22											
*05:44		2											13																				
*05:45, 05:169, 05:182 ⁴	1	2										11																					
*05:46	1		3														17												29				
*05:47	1		3													16													29				
*05:49	1	2											13																				
*05:51	1	2																												30			
*05:52	1				5		7																						29				
*05:54		2																															
*05:56	1			4																								28					
*05:57	1		3																		20								29				
*05:60				4																19													
*05:62	1	2																			20												
*05:63	1	2																				21											
*05:66:01-05:66:02	1			4		6																											
*05:69	1	2												14																			
*05:71	1	2								10														24									
*05:72		2															18																
*05:73		2											14										24										
*05:75	1	2																							25								
*05:78	1			4																							26						
*05:79	1		3																									27	29				
*05:80	1	2															18					24											
*05:82																						21											
*05:83, 05:175			3																										29				
*05:87Q	1		3					8																					29				
*05:88	1	2																				21						27					
*05:90N	1		3									12					17												29				
*05:96	1			4							11	12																					
*05:97	1																							23									
*05:98				4										14										24									
*05:102	1		3																									29	31				
*05:103	1	2																									26						
*05:104	1	2													15							21											
*05:105																16														30			
*05:106	1		3														17				20								29				
*05:107	1	2																										28					
*05:108	1			4													17																
*05:110N, 05:156 ⁴	1	2																										27					
*05:111	1	2																								25			30				
*05:112	1	2					7																		25								
*05:116			3											14										24					29				
*05:128N		2											13				18																
*05:130				4									13																				
*05:132Q	1				5		8																						29				
*05:134	1																																
*05:143	1			4																				23									
*05:145, 05:153	1		3											14															29				
*05:218		2		4																													
*05:231	1		3																								26		29				
*05:236N	1			4																								27					
*05:240	1	2		4		6																											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*02:64, 02:79, 03:21, 06:243, 06:255												11																					
*02:96N												12					17																
*02:99, 03:249, 03:312, 06:105, 06:185, 06:225																16																	

Negative Control

