

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
DQB1 allele ^{1,2,3}																																	
*05:01:01:01-05:01:08, 05:01:11-05:01:41, 05:27 , 05:48, 05:55, 05:61, 05:68, 05:70, 05:74, 05:76, 05:81, 05:84, 05:89:01-05:89:02, 05:92-05:93, 05:95, 05:99, 05:114-05:115, 05:120, 05:122, 05:124, 05:126-05:127, 05:133, 05:137-05:139, 05:141, 05:144, 05:148, 05:150-05:152, 05:154-05:155, 05:157-05:158, 05:160, 05:162-05:164, 05:166-05:168, 05:171, 05:173, 05:176-05:177, 05:180, 05:183-05:185N, 05:187-05:188, 05:190, 05:193-05:195, 05:197, 05:215N-05:217, 05:219, 05:223, 05:225-05:226, 05:228, 05:230, 05:232, 05:234, 05:237, 05:246, 05:248-05:249, 05:252	1	2																															
*05:01:09-05:01:10	1	2		4																													
*05:02:01:01-05:02:01:07, 05:02:03-05:02:12, 05:02:14-05:02:17, 05:02:19-05:02:23, 05:53, 05:59, 05:65, 05:86, 05:94, 05:100, 05:113, 05:118, 05:123, 05:125, 05:135-05:136, 05:142, 05:165, 05:170, 05:172, 05:174, 05:178-05:179, 05:181, 05:186, 05:189, 05:192, 05:196, 05:198-05:199, 05:222, 05:227, 05:229, 05:238-05:239, 05:241, 05:247, 05:250-05:251, 05:253-05:254	1		3																											29			
*05:02:02, 05:02:18	1		3	4																										29			
*05:02:13	1		3																														
*05:03:01:01-05:03:05, 05:03:07-05:03:13, 05:03:15-05:03:19, 05:03:21-05:03:23, 05:43:02, 05:64, 05:85, 05:91, 05:109, 05:121, 05:129, 05:140, 05:147, 05:149, 05:161, 05:191, 05:200, 05:202-05:214, 05:220-05:221, 05:224N, 05:233, 05:235N, 05:243, 05:245	1			4																													
*05:03:06	1			4											15																		
*05:03:14, 05:03:20	1														15						22												
*05:04, 05:77, 05:146	1				5																									29			
*05:05:01-05:05:02	1		3			6																								29			
*05:06:01-05:06:02, 05:50	1			4			7																										
*05:07	1	2					7																										
*05:08	1			4					9																								
*05:09	1			4				8																									
*05:10	1			4	5																												
*05:11:01-05:11:02	1	2				6																											
*05:12, 05:242 ⁴	1	2								10																							
*05:13, 05:201	1			4							11																						
*05:14	1		3									12																		29			
*05:15	1			4									13																				
*05:16, 05:101 ⁴	1			4										14																			
*05:17	1		3												15															29			
*05:18	1	2							9																								
*05:19, 05:97 ⁵	1		3																			23							29				
*05:20, 05:159 ⁴	1	2													16																		
*05:21																		19															
*05:22	1	2																				23											
*05:23	1									10																							
*05:24	1			4														18															
*05:25	1	2										12				17					21												
*05:26	1		3																										29				
*05:28	1			4															19														
*05:29	1	2						8																									
*05:30	1	2																	19														

Negative Control

[illegible]

[illegible]