

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*05:19	1		3																			23							29				
*05:20:01, 05:159 ⁴	1	2														16																	
*05:21																			19														
*05:22	1	2																					23										
*05:23	1									10																							
*05:24	1			4														18															
*05:25	1	2										12					17																
*05:26	1		3																			21								29			
*05:27, 05:29 ⁴	1	2						8																									
*05:28	1			4																19													
*05:30	1	2																	19														
*05:31	1	2															17																
*05:32	1	2									11										20												
*05:33	1		3										13																		29		
*05:34			3										13																		29		
*05:35, 05:117 ⁴	1		3								11																				29		
*05:36	1		3					9																							29		
*05:37	1		3												15							21									29		
*05:38, 05:119	1			4																	20												
*05:39	1			4																		21											
*05:40	1			4									13	14																			
*05:41N	1			4								12					17																
*05:42	1			4								11										20											
*05:43:01, 05:131, 05:276	1																						22										
*05:44		2												13																			
*05:45, 05:169, 05:182 ⁴	1	2										11																					
*05:46	1		3														17														29		
*05:47	1		3													16															29		
*05:49	1	2												13																			
*05:51	1	2																													30		
*05:52	1				5		7																								29		
*05:54		2																													31		
*05:56	1			4																										28			
*05:57	1		3																			20									29		
*05:58	1			4										14					18														
*05:60				4																19													
*05:62	1	2																				20											
*05:63	1	2																					21										
*05:66:01-05:66:02	1			4		6																											
*05:67	1			4																											31		
*05:69	1	2													14																		
*05:71	1	2									10														24								
*05:72		2																	18														
*05:73		2												14										24									
*05:75	1	2																								25							
*05:78	1			4																										26			
*05:79	1		3																										27		29		
*05:80	1	2																18						24									
*05:82																							21										
*05:83, 05:175			3																												29		
*05:84	1	2										12																					
*05:87Q	1		3					8																							29		
*05:88	1	2																					21							27			
*05:90N	1		3									12					17														29		
*05:96	1			4								11	12																				
*05:97	1																							23									
*05:98				4										14										24									
*05:102	1		3																											29	31		
*05:103	1	2																												26			
*05:104	1	2													15								21										
*05:105																16																30	
*05:106	1		3														17					20									29		
*05:107	1	2																												28			
*05:108	1			4													17																
*05:110N, 05:156 ⁴	1	2																												27			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*05:111	1	2																														30	
*05:112	1	2					7																										

Negative Control

*05:116		3				14		24		29	
*05:128N	2					13	18				
*05:130		4				13					
*05:132Q	1		5	8						29	
*05:143	1	4						23			
*05:145, 05:153	1	3				14				29	
*05:211	1		4			12					
*05:218		2	4								
*05:231	1	3							26	29	
*05:236N	1		4							27	
*05:240	1	2	4	6							
*02:64, 02:79, 03:21, 03:452, 06:243, 06:255						11					
*02:96N						12		17			
*02:99, 03:249, 03:312, 06:105, 06:185, 06:225, 06:407							16				
*02:157, 03:190, 04:42, 06:223, 06:406											31
*02:169, 03:319, 03:347, 04:70, 06:103, 06:284, 06:350						12					
*03:03:05							15				30
*03:03:28, 03:10:02:01-03:10:03, 03:38:02, 03:183, 06:01:12, 06:02:43, 06:03:33, 06:51:01, 06:66, 06:96:01, 06:99:01, 06:118:04, 06:172, 06:322:01								22			
*03:226, 03:272, 04:28		4									
*03:309			5								
*03:317:01-03:317:02, 06:111, 06:265								17			
*03:420				8							
*03:432								20			
*04:02:07									w		
*04:56						11		20			
*06:02:19, 06:03:08, 06:09:11, 06:13:03											30
*06:06			?								
*06:28, 06:56, 06:79:01, 06:89									24		
*06:79:02									24		30
*06:129, 06:241			6								
*06:146:01			6					20			
*06:156, 06:162, 06:169		4					14				
*06:158N										27	
*06:221, 06:273, 06:326				9							

Negative Control