

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
DQB1 allele ^{1,2,3}																																
*05:01:01:01-05:01:08, 05:01:11-05:01:45, 05:48, 05:55, 05:61, 05:68, 05:70, 05:74, 05:76, 05:81, 05:89:01-05:89:02, 05:92-05:93, 05:95, 05:99, 05:114-05:115, 05:120, 05:122, 05:124, 05:126-05:127, 05:133, 05:137-05:139, 05:141, 05:144, 05:148, 05:150-05:152, 05:154-05:155, 05:157-05:158, 05:160, 05:162-05:164, 05:166-05:168, 05:171, 05:173, 05:176-05:177, 05:180, 05:183-05:185N, 05:187-05:188, 05:190, 05:193-05:195, 05:197, 05:215N-05:217, 05:219, 05:223, 05:225-05:226, 05:228, 05:230, 05:232, 05:234, 05:237, 05:246, 05:248-05:249, 05:252, 05:255-05:256, 05:258, 05:261, 05:263, 05:266-05:272, 05:275, 05:277, 05:279, 05:283N, 05:285-05:286, 05:288, 05:290-05:291	1	2																														
*05:01:09-05:01:10, 05:292	1	2		4																												
*05:02:01:01-05:02:01:12, 05:02:03-05:02:12, 05:02:14-05:02:17, 05:02:19-05:02:26, 05:53, 05:59, 05:65, 05:86, 05:94, 05:100, 05:113, 05:118, 05:123, 05:125, 05:135-05:136, 05:142, 05:165, 05:170, 05:172, 05:174, 05:178-05:179, 05:181, 05:186, 05:189, 05:192, 05:196, 05:198-05:199, 05:222, 05:227, 05:229, 05:238-05:239, 05:241, 05:247, 05:250-05:251, 05:253-05:254, 05:257, 05:274, 05:278, 05:284, 05:287	1		3																										29			
*05:02:02, 05:02:18	1		3	4																									29			
*05:02:13	1		3																													
*05:03:01:01-05:03:05, 05:03:07-05:03:13, 05:03:15-05:03:19, 05:03:21-05:03:23, 05:03:25-05:03:28, 05:43:02, 05:64, 05:85, 05:91, 05:109, 05:121, 05:129, 05:140, 05:147, 05:149, 05:161, 05:191, 05:200, 05:202-05:210, 05:212-05:214, 05:220-05:221, 05:224N, 05:233, 05:235N, 05:243, 05:245, 05:260, 05:264-05:265N, 05:273N, 05:280-05:282, 05:289, 05:293-05:294	1			4																												
*05:03:06	1			4											15																	
*05:03:14, 05:03:20	1														15						22											
*05:03:24, 05:134, 05:244	1																															
*05:04, 05:77, 05:146, 05:262	1				5																								29			
*05:05:01-05:05:02	1		3			6																							29			
*05:06:01-05:06:02, 05:50	1			4			7																									
*05:07, 05:20:02	1	2					7																									
*05:08	1			4					9																							
*05:09	1			4				8																								
*05:10	1			4	5																											
*05:11:01-05:11:02	1	2				6																										
*05:12, 05:242 ⁴	1	2								10																						
*05:13, 05:201, 05:259 ⁴	1			4							11																					
*05:14	1		3									12																	29			
*05:15	1			4									13																			
*05:16, 05:101 ⁴	1			4										14																		
*05:17	1		3												15														29			
*05:18:01-05:18:02	1	2							9																							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32

Negative Control

[illegible]

*05:130			4					13									
*05:132Q	1			5		8										29	
*05:143	1		4									23					
*05:145, 05:153	1	3						14								29	
*05:211	1		4					12									
*05:218		2	4														
*05:231	1		3											26		29	
*05:236N	1		4											27			
*05:240	1	2	4		6												
*02:64, 02:79, 03:21, 03:452, 06:243, 06:255								11									
*02:96N								12			17						
*02:99, 03:249, 03:312, 06:105, 06:185, 06:225, 06:407										16							
*02:157, 03:190, 04:42, 06:223, 06:406																	31
*02:169, 03:319, 03:347, 04:70, 06:103, 06:284, 06:350								12									
*03:03:05										15							30
*03:03:28, 03:10:02:01-03:10:03, 03:38:02, 03:183, 06:01:12, 06:02:43, 06:03:33, 06:51:01, 06:66, 06:96:01, 06:99:01, 06:118:04, 06:172, 06:322:01												22					
*03:226, 03:272, 04:28			4														
*03:309				5													
*03:317:01-03:317:02, 06:111, 06:265											17						
*03:420						8											
*03:432											20						
*04:02:07												w					
*04:56								11			20						
*06:02:19, 06:03:08, 06:09:11, 06:13:03																	30
*06:06					?												
*06:28, 06:56, 06:79:01, 06:89													24				
*06:79:02													24				30
*06:129, 06:241					6												
*06:146:01					6						20						
*06:156, 06:162, 06:169			4						14								
*06:158N															27		
*06:221, 06:273, 06:326						9											

Negative Control