

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		
DQA1 allele <sup>1,2</sup>																																		
*01:01:01-01:01:02	1	2						8					13											24										
*01:02:01:01-01:02:03		2	3					8					13											24										
*01:02:04			2	3									13											24										
*01:03:01:01-01:03:01:02				3	4			8					13											24										
*01:04:01:01-01:04:02	1	2				6							13											24										
*01:05	1	2			5								13											24										
*01:06		2	3			?	7	?					13											24										
*01:07	1	2				6							13											24	25									
*01:08		2	3					8					13											24	25					29				
*01:09		2	3					8					13									21		24										
*01:10			3	4				8					13									22		24										
*01:11		2	3					8					13											24			27							
*01:12	1	2						8					13											24						29				
*02:01								8	9				13											w										
*03:01:01								8		10			13											24										
*03:02								8		10	11	12												24										
*03:03:01-03:03:02								8		10		12												24										
*04:01:01								8					13	14										23	24									
*04:01:02:01-04:01:02:02								8					14											23	24									
*04:02								8					13	14										23	24		26							
*04:03N								8					13											23	24			27						
*04:04								8					13	14										23	24				28					
*05:01:01:01-05:01:02								8					13	14	15		17		19					24										
*05:02								8					13	14	15	16	?	?	19	?				?										
*05:03								8					13	14	15				18	19				24										
*05:04								8					13	14	15		?	?		20	?			?										
*05:05:01:01-05:05:01:03								8					13	14	15		17		19		21													
*05:06								8					13	14	15				18	19				24							30			
*05:07								8					13	14	15	16			18	19				24										
*05:08								8					13	14	15		17		19		21											31		
*05:09								8					13	14	15		17		19		21					25								
*05:10								8					13	14					17		?			?		26								
*05:11								8					13	14	15		17		19		21											30		
*06:01:01-06:01:02								8					13										22	23	24									
*06:02								8					13										22	23	24						29			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32		

Negative Control