

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
DQA1 allele ^{1,2,3}																																	
*01:01:01-01:01:03	1	2						8					13											24									
*01:02:01:01-01:02:03		2	3					8					13											24									
*01:02:04			2	3									13											24									
*01:03:01:01-01:03:01:02				3	4			8					13											24									
*01:04:01:01-01:04:02	1	2			5	6							13											24									
*01:05:01-01:05:02	1	2			5								13											24									
*01:06		2	3		?	?	7	?					13											24									
*01:07Q	1	2			5	6							13											24	25								
*01:08		2	3		?	?		?					13											24					29				
*01:09		2	3		?	?		?					13								21			24									
*01:10			3	4				8					13									22		24									
*01:11		2	3					8					13											24			27						
*01:12	1	2			?	?		?					13											24					29				
*01:13		2	3					8					13											24	25								
*02:01								8	9				13											w									
*03:01:01, 03:01:03								8		10			13											24									
*03:02								8		10	11	12												24									
*03:03:01-03:03:02								8		10		12												24									
*04:01:01								8					13	14										23	24								
*04:01:02:01-04:01:02:02								8						14										23	24								
*04:02								8					13	14										23	24	26							
*04:03N								8					13											23	24		27						
*04:04								8					13	14										23	24				28				
*05:01:01:01-05:01:02								8					13	14	15		17		19					24									
*05:02								8					13	14	15	16	?	?	19		?			?									
*05:03								8					13	14	15				18	19				24									
*05:04								8					13	14	15		?	?		20	?			?									
*05:05:01:01-05:05:01:03								8					13	14	15		17		19		21												
*05:06								8					13	14	15				18	19				24						30			
*05:07								8					13	14	15	16			18	19				24									
*05:08								8					13	14	15		17		19		21										31		
*05:09								8					13	14	15		17		19		21					25							
*05:10								8					13	14			17				?		?			26							
*05:11								8					13	14	15		17		19		21									30			
*06:01:01-06:01:02								8					13										22	23	24								
*06:02								8					13										22	23	24					29			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	

Negative Control