

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
DQA1 allele^{1,2,3}																																			
*01:01:01:01-01:01:01:03, 01:01:01:05-01:01:06, 01:01:08-01:01:11, 01:27, 01:49	1	2						8					13												24								31		
*01:01:07	1	2						8					13												24										
*01:02:01:01-01:02:03, 01:02:05-01:02:10, 01:19, 01:21, 01:23, 01:25, 01:38:01:01-01:38:01:02, 01:41, 01:51, 01:54, 01:58		2	3					8					13												24										
*01:02:04		2	3										13												24										
*01:03:01:01-01:03:01:12, 01:03:03-01:03:05, 01:30, 01:44-01:45, 01:57			3	4				8					13												24										
*01:03:02			3	4	?			8					13												24										
*01:04:01:01-01:04:02, 01:53, 01:64	1	2			5	6							13												24										
*01:04:03	1	2			5	6		8					13												24										
*01:05:01:01-01:05:02, 01:22	1	2			5								13												24								31		
*01:05:03		2			5								13												24								31		
*01:06		2	3		?	?	7	?					13												24								?		
*01:07Q	1	2			5	6							13												24	25									
*01:08		2	3		?	?		?					13												24		26						?		
*01:09		2	3		?	?		?					13									21			24								?		
*01:10			3	4				8					13		15										24										
*01:11		2	3					8					13												24										
*01:12	1	2			5								13												24		26						31		
*01:13		2	3		?	?		?					13												24	25							?		
*01:14			3	4				8					13												24				28						
*01:15N			3	4	?	?		?					13						19						24								?		
*01:16N		2	3					8					13												24		26								
*01:17, 01:33			3	4	?			?					13												24										
*01:18, 01:26, 01:34-01:35	1	2			?	?		?					13												24								?		
*01:20, 01:28, 01:31-01:32, 01:40Q		2	3		?			?					13												24										
*01:24			3	4				8					13												24	25									
*01:29	1	2			?	6		?					13												24										
*01:36		2	3		?	?		?					13												24									?	
*01:37	1	2			?			?					13												24									31	
*01:39		2	3		5								13												24										
*01:42		2	3		?	?		?					13		15										24									?	
*01:43	1	2						8					13												24	25								31	
*01:46		2	3		5	?		8					13												24									?	
*01:47			3	4	5	?		8					13												24									?	
*01:48		2	3					8																	24										
*01:50			3	4	?	?		8					13												24									?	
*01:52, 01:62		2	3		?	?		8					13												24									?	
*01:55, 01:59	1	2			?			8					13												24									31	
*01:56, 01:60-01:61	1	2			?	?		8					13												24									?	
*01:63		2	3		?			8					13												24									31	
*01:65			3	4	5								13												24										
*02:01:01:01								8	9				13												w										
*02:01:01:02, 02:01:03- 02:01:06, 02:03-02:05, 02:07-02:15, 02:18								8	9				13													24									
*02:01:02								8	9																24										
*02:02N								8	9				13										22		24										
*02:16								8					13												24										
*02:17								8	9			12													24										
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32			
*03:01:01:01, 03:01:03- 03:01:06, 03:06, 03:10, 03:12, 03:15								8		10			13												24										

Negative Control

