

DQB1*03:18, 03:37, 03:63, 03:66N-03:68, 03:70, 03:106-03:107, 03:125, 03:146, 03:153, 03:161, 03:174-03:175, 03:178-03:179, 03:185, 03:203-03:204, 03:210, 03:213N, 03:215, 03:221, 03:223-03:224, 03:228-03:229, 03:233, 03:237N, 03:261, 03:269N, 03:273-03:274, 03:278, 03:352														14	16	17	19	?	22				
DQB1*03:19:01:01-03:19:04, 03:150, 03:286, 03:376N, 03:394, 03:401														14			19	20	21				
DQB1*03:24, 03:164, 03:329, 03:347																	19	21					
DQB1*03:48, 03:52-03:54, 03:56-03:60, 03:69, 03:71, 03:75-03:77, 03:82, 03:101-03:103, 03:108, 03:118N, 03:120-03:122, 03:127-03:131, 03:133-03:135, 03:140, 03:142-03:144, 03:147-03:148, 03:152, 03:157-03:160, 03:162-03:163, 03:166-03:167, 03:170-03:173, 03:183, 03:187-03:188, 03:192-03:194, 03:201-03:202, 03:207-03:208, 03:216, 03:218-03:219, 03:232, 03:235-03:236, 03:255-03:257, 03:260, 03:267-03:268, 03:275, 03:353, 03:387														14			19	?	21				
DQB1*03:62, 03:81, 03:225														14	16	17	19	?					
DQB1*03:74, 03:111, 03:141, 03:356N														14	15		19	?					
DQB1*03:79, 03:97, 03:319															15	17	19						
DQB1*03:104, 03:258														14	15	17		?					
DQB1*03:109, 03:154														14				?					
DQB1*03:165																	19	?	21				
DQB1*03:189														14	16	17		?	22				
DQB1*03:308, 03:413														14	16	17	19						
DQB1*03:335														14				20					
DQB1*03:349															16	17	19		22				
DQB1*03:381																	19	20	21				
DQB1*04:01:01:01-04:01:02, 04:01:04-04:02:01:01, 04:02:01:04-04:02:15, 04:02:17-04:02:20, 04:04-04:08, 04:10-04:43, 04:45-04:62, 04:64-04:69, 04:72-04:74, 04:76-04:77, DQB1*05:01:01:01-05:02:10, 05:02:12-05:13, 05:15-05:83, 05:85-05:195, 05:197-05:210, 05:212-05:245, DQB1*06:01:01:01-06:02:04, 06:02:06-06:02:42, 06:03:01:01-06:03:09, 06:03:11-06:03:25, 06:03:27-06:03:32, 06:03:35-06:04:06, 06:04:08-06:18:02, 06:20-06:28, 06:30-06:37, 06:39-06:50, 06:51:02-06:62, 06:64-06:65, 06:67-06:85, 06:88-06:95, 06:96:02-06:101, 06:105-06:122, 06:124-06:138, 06:140-06:167, 06:169-06:171, 06:173-06:189, 06:190:02-06:217, 06:219, 06:221-06:226, 06:228-06:245, 06:247-06:283, 06:285-06:319, 06:321, 06:323-06:336, 06:338-06:349, 06:351-06:362														14									
DQB1*04:03:01-04:03:02, DQB1*06:02:05, 06:03:26, 06:03:34, 06:04:07, 06:19:01-06:19:02, 06:87, 06:190:01														14		17							
DQB1*06:02:43, 06:03:33, 06:51:01, 06:66, 06:96:01, 06:172														14	15								
DQB1*06:29, 06:123														14	16				22				
DQB1*06:63														14	16	17							
DQB1*06:139, 06:337														14	16	17			22				
DQB1*06:246														14			18						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23

Negative Control