

DQB1 high 2019 4H0 R04

DQB1 high 2019 4HO R04

DQB1 high 2019 4H0 R04

*06:184	+	?			?	+++		?	+		?						+	+	
*06:186	+	w				++	w												
*06:187	+				++	++		+	+								+		
*06:188	+				++	++		+	+										
*06:190:01-06:190:02	+	?			?	++		?	+		?						++		
*06:191	+	?			?	++		?	+		?							++	
*06:193N	+	?			?	++		?	+		?								
*06:194	+	?		+	?			+	?		?								
*06:195	+	?			?	++		?	+		?								
*06:197	+	?			?	++		?			?								
*06:202	+	?			?		++	?			?								
*06:206:01	+	?			?		+	++	?	+									
*06:206:02	+	w					+	++w	+	+									
*06:209	+			+		+			+								+	+	
*06:210	+	?		+	?	++		?			?								
*06:212	+	?			?		+	w	?	+									
*06:214	+	?		+	?			+	?		?								
*06:216N	+				++		++		w		++								
*06:217	+	w																	
*06:218	+				++														
*06:221	+				++														
*06:224	+				++														
*06:227					++														
*06:228	+				++														
*06:231	+	?			?		++		?										+
*06:233	+	?		+	?		?		?										
*06:238	+	?			++		?		+										
*06:239	+			+															
*06:241	+	+					++		?				?	+					
*06:243	+				+				+										
*06:244	+	?				+	+		?				?						
*06:245, 06:251	+	?			+				?										
*06:246	+	?		+	+				?				?						
*06:247	+	?		+				+	?				+	?					
*06:248	+	?					++		?				+						
*06:249	+	?			++		?		+				?						
*06:250	+				++				+										
*06:252N	+	?				++		?				+	?						+
*06:254	+	w				++		w	w				+						
*06:255	+				++			+											
*06:260	+	?		+			+		?			?							
*06:263, 06:268, 06:277	+	?		+			+		?			?							
*06:265	+	w				++		w	w				+						
*06:269	+	?			++		+	?				+							
*06:270	+				++							+							
*06:273	+				++			+											
*06:275	+	w			++			w					+						
*06:276	+				++			+											
*06:278	+				++			+											
*06:280, 06:283	+	?				++		?					+						
*06:281-06:282	+	?				++		?					+						
*06:284					++			+											
Well No.	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9
DQB1 allele ^{1,2,3}	1	2	3	4	5	6	7	8	9	0	1	2	3	4	5	6	7	8	9