

Well No.	1	2	3	4	5	6	7	8	
DQA1, DQB1 or DRB1 alleles^{1,2}									
DQA1*01:01:01-01:01:01:02, 01:04:01:01-01:05:02, 01:07, 01:12	1	2						Negative Control	
DQA1*01:02:01-01-01-02:04 , 01:06, 01:08-01:09, 01:11		2	3						
DQA1*01:03:01:01-01:03:01:02, 01:10			3	4					
DQB1*03:30, DQB1*03:72, DQB1*03:100, DQB1*03:132, DQB1*04:01:01-04:02:04, DQB1*04:02:06, DQB1*04:02:08-04:03:02, DQB1*04:06-04:08, DQB1*04:12-04:14, DQB1*04:16-04:27, DQB1*05:38, DQB1*06:139						6			
DQB1*04:10, DQB1*06:02:05 , DQB1*06:19:01					5	6			
DQB1*06:02:01-06:02:04, DQB1*06:02:06-06:02:23 , DQB1*06:14:01-06:16, DQB1*06:19:02-06:20, DQB1*06:23-06:24, DQB1*06:33, DQB1*06:37, DQB1*06:46-06:50, DQB1*06:51:02, DQB1*06:68, DQB1*06:70-06:84, DQB1*06:95, DQB1*06:97, DQB1*06:107, DQB1*06:109, DQB1*06:111-06:117, DQB1*06:122, DQB1*06:124-06:127, DQB1*06:136-06:138, DQB1*06:146-06:147, DQB1*06:150-06:152, DQB1*06:156, DQB1*06:159, DQB1*06:161-06:163, DQB1*06:166					5				
DRB1*15:01:01-01-15:114, DRB1*15:116-15:117							7		
Well No.	1	2	3	4	5	6	7		8