

Well No.	1	2	3	4	5	6	7	8
<b>DQA1, DQB1 or DRB1 alleles<sup>1,2</sup></b>								
DQA1*01:01:01:01-01:01:01:03, DQA1*01:04:01:01-01:05:02, DQA1*01:07Q, DQA1*01:12	1	2						
<b>DQA1*01:02:01:01-01:02:04</b> , 01:06, 01:08-01:09, 01:11, 01:13		2	3					
DQA1*01:03:01:01-01:03:01:06, 01:10			3	4				
DQB1*03:30, DQB1*03:72, DQB1*03:100, DQB1*03:132, DQB1*03:215, DQB1*04:01:01-04:02:04, DQB1*04:02:06, DQB1*04:02:08-04:03:02, DQB1*04:06-04:08, DQB1*04:12-04:14, DQB1*04:16-04:34, DQB1*05:38, DQB1*05:119, DQB1*06:139						6		
DQB1*04:10, <b>DQB1*06:02:05</b> , DQB1*06:19:01					5	6		
<b>DQB1*06:02:01-06:02:04, DQB1*06:02:06-06:02:25</b> , DQB1*06:14:01-06:16, DQB1*06:19:02-06:20, DQB1*06:23-06:24, DQB1*06:33, DQB1*06:37, DQB1*06:46-06:50, DQB1*06:51:02, DQB1*06:68, DQB1*06:70-06:84, DQB1*06:95, DQB1*06:97, DQB1*06:107, DQB1*06:109, DQB1*06:111-06:117, DQB1*06:122, DQB1*06:124-06:127, DQB1*06:136-06:138, DQB1*06:146:01-06:147, DQB1*06:150-06:152, DQB1*06:156, DQB1*06:159, DQB1*06:161-06:163, DQB1*06:166, DQB1*06:173-06:175, DQB1*06:178-06:179N, DQB1*06:182-06:183, DQB1*06:188, DQB1*06:192, DQB1*06:197-06:198, DQB1*06:200-06:201					5			
<b>DRB1*15:01:01-15:131</b>							7	
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>

Negative Control