

Well No.	DR	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
DRB1 allele ¹	ser ²																				
*010101-010204, 0104-0119	1	1																			
*0103	103	2																			
*030101-030106, 0304-0306, 0309, 0311-0316, 0318-0320, 0322, 0323, 0325, 0326, 0328, 0330, 0331, 0333, 0334, 0336, 0337	17					5	6											17			
*030201-0303, 0327, 0329, 0338	18					5	7											17			
*0307, 0317, 0321, 0324, 0332, 0335, 0339	-					5												17			
*0308	-					5	6				11			13						19	
*0310	-					5	6													19	20
*040101-0461, 0463-0468, 0470-0472	4								8												
*0462, 0469, 0473	-								8										18		
*070101-070102, 0703-0715	7									9											
*080101-080203, 080401-0807, 0811, 0816, 0817, 0822, 0824, 0826, 0828	8										10					16					
*080302, 0810, 0812-0815, 0818, 0819, 0823, 0825, 0827, 0829, 0830, 0833, 0834	8										10										
*0808	-										10					16			18		
*0809, 0821, 1415	8										10					16					20
*0820, 1318, 1347, 1355	8,13														15	16		17			20
*0831, 1167	-										10			13		16					
*0832	-										10				14						20
*090102-0902, 0904-0906	9											11									19
*0903	9											11									
*100101-1002	10												12								
*110101-110107, 110401-110403, 110601-110602, 111201, 1115, 1124, 112701-1130, 1132, 1133, 1135, 1137-1139, 1143, 1144, 1147, 1149-1151, 115401-115402, 1156, 1160-1162, 1166	11													13		15	16				
*110201-1103, 111101-111102, 111401-111402, 1121, 1136, 1141, 1148, 1163, 116501-116502, 1168	11						6							13		15	16				
*110404, 1109, 1110, 11202, 1146, 1158	11													13		15	16			19	
*1105	11													13		16					
*1107, 1153	-					5					11			13							
*110801-110802, 1118-111902, 1142, 1157	11													13		15					
*1113, 1117, 1152	11													13						19	20
*1116, 1120, 1140, 1159	-						6							13		15	16			19	
*1122	-													13						19	
*1123, 1125	11													13		15	16				20
*1126, 1134	11							7						13							
*1131, 1145, 1164	-													13		15					20
*1155	-													13		16					20
*120101-120102, 120302, 1205-1212, 1214	12														14						
*120201-120202, 1213, 1215	12														14		16				
*1204	-										10				14					19	



