

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
<b>DRB1 allele<sup>1,2,3</sup></b>																								
*01:01:01-01:01:29, 01:38, 01:41, 01:45, 01:47-01:49, 01:53, 01:56, 01:59-01:60, 01:64-01:66, 01:69-01:73, 01:75, 01:78	1		3																					
*01:02:01-01:02:05, 01:02:07-01:02:12, 01:43, 01:46, 01:57, 01:76			2	3																			23	
*01:02:06, 01:58				3																				23
*01:03	1			4			7																	
*01:04			3		5				9															
*01:05	1		3			6																		
*01:06			3						9	10														
*01:07	1		3					8																
*01:08	1		3						9															
*01:09	1		3							10														
*01:10	1		3	4								11												
*01:11:01-01:11:02	1		3		5																			
*01:12, 01:37	1		3									12												
*01:13	1		3										13											
*01:14	1		3											14										
*01:15	1						7			10														
*01:16	1		3		5										15									
*01:17	1		3				7																	
*01:18	1		3								11													
*01:19	1		3													16								
*01:20:01			3						9															
*01:20:02			2	3					9															
*01:21	1		3												15	16								
*01:22	w		w														17							
*01:23		2											13	14										23
*01:24:01-01:24:02	1													14		16								
*01:25	1		3															18						
*01:26		2	3								12													23
*01:27	1		3					8																
*01:28	1		3										13		15									
*01:29:01	1		3		5	6																		
*01:29:02	1					6																		
*01:30	1		3																					
*01:31	1		3			6					12													
*01:32	1		3																				21	
*01:33N																					20			
*01:34		2	3																			21	23	
*01:35			3		5																		23	
*01:36	1		3																				22	
*01:39N	1			4			7																	19
*01:40N, 01:67 <sup>4</sup>	1		3																					19
*01:42				4			7																	23
*01:44:01-01:44:02		2	3								11													23
*01:50	1		3					8																18
*01:51	1												13	14										
*01:52N	1		3							10														22
*01:54		2												14		16								23
*01:55	1				5										15									
*01:61				4			7		9															
*01:62N	1		3	4																			21	
*01:63			3																					
*01:68N	1		3								11			14										
*01:74		2	3		5	6																		23
*01:77	1		3																					19
*13:165																								17
*14:112								8																
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>

Negative Control