

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
<b>DRB1 allele<sup>1,2,3</sup></b>																								
*01:01:01:01-01:01:35, 01:01:37-01:01:42, 01:01:44-01:01:46, 01:38, 01:41, 01:45, 01:47-01:49, 01:53, 01:56, 01:59-01:60, 01:62:02N-01:66, 01:69-01:73, 01:75, 01:78-01:80, 01:82, 01:84, 01:91Q-01:92, 01:95, 01:98, 01:100, 01:103-01:110, 01:113, 01:115, 01:117, 01:120-01:121, 01:124-01:125, 01:128-01:130, 01:132-01:133, 01:135	1		3																					
*01:01:36	1	2	3																					
*01:01:43	1		3										13	14										
*01:02:01:01-01:02:05, 01:02:07-01:02:13, 01:43, 01:46, 01:57, 01:76, 01:81, 01:83, 01:86-01:88, 01:93-01:94, 01:96-01:97, 01:122-01:123, 01:126, 01:136			2	3																				23
*01:02:06, 01:58, 01:116				3																				23
*01:03:01:01-01:03:05, 01:102, 01:127	1			4			7																	
*01:04			3		5				9															
*01:05	1		3			6																		
*01:06			3						9	10														
*01:07	1		3					8																
*01:08	1		3						9															
*01:09	1		3							10														
*01:10	1		3	4							11													
*01:11:01-01:11:02	1		3		5																			
*01:12, 01:37	1		3									12												
*01:13	1		3										13											
*01:14	1		3											14										
*01:15	1						7			10														
*01:16	1		3		5										15									
*01:17	1		3				7																	
*01:18:01-01:18:02	1		3								11													
*01:19	1		3																					
*01:20:01, 01:85			3						9															
*01:20:02			2	3					9															
*01:21	1		3												15	16								
*01:22	1		w																					
*01:23, 01:111		2											13	14										23
*01:24:01-01:24:02	1													14		16								
*01:25	1		3																					
*01:26		2	3									12												23
*01:27	1		3					8																
*01:28	1		3										13	15										
*01:29:01	1		3		5	6																		
*01:29:02	1					6																		
*01:30	1		3																					
*01:31	1		3				6					12												
*01:32	1		3																					21
*01:33N	1																							20
*01:34		2	3																					21
*01:35			3		5																			23
*01:36	1		3																					22
*01:39N	1		4				7																	19
*01:40N, 01:67 <sup>d</sup>	1		3																					19
*01:42				4			7																	23
*01:44:01-01:44:02		2	3							11														23
*01:50	1		3					8																18
*01:51	1												13	14										
*01:52N	1		3							10														22
*01:54		2												14	16									23
*01:55	1				5										15									
*01:61				4			7		9															
*01:62:01N	1		3	4																				21
*01:68N	1		3							11			14											
*01:74		2	3		5	6																		23
*01:77, 01:119	1		3																					19
*01:89	1			4	5	7																		
*01:90	1														15									
*01:99	1			4																				
*01:101			3		5				9						15									
*01:112, 01:131N, 01:134	1																							
*01:114		2	3				7																	23
*01:118		2	3		5	6	7																	23
*10:40, DRB5*01:60																								18
*13:165, DRB3*02:63																								17
*14:112								8																
*14:196																								19
DRB3*01:47								8																17
DRB3*03:01:08			3		5																			
DRB4*01:132																								19
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>

Negative Control