

INTERPRETATION TABLE				
DRB1*10 SSP subtyping				
Amplification patterns of the DRB1*10 alleles				
	Well			
	1	2	3	4
Length of spec.	210	205	205	70
PCR product				
Length of int.	515	515	430	430
pos. control¹				
5'-primer²	31	31	31	31
	5'-gC g ^{3'}	5'-gC g ^{3'}	5'-gC g ^{3'}	5'-gC g ^{3'}
3'-primer³	87	86	86	40
	5'-CTC ^{3'}	5'-C AC ^{3'}	5'-C CA ^{3'}	5'-C gA ^{3'}
Well No.	1	2	3	3
DRB1 allele				
*100101-100103	1	2		
*1002	1		3	
*1003	1	2		4
DRB1 allele				
Well No.	1	2	3	3