

INTERPRETATION TABLE					
DRB1*10 SSP subtyping					
Amplification patterns of the DRB1*10 alleles					
	Well				
	1	2	3	4	5
Length of spec.	210	205	205	70	160
PCR product					
Length of int.	515	515	430	430	430
pos. control¹					
5'-primer²	31(178)	31(178)	31(178)	31(178)	31(178)
	5' -gCg 3'	5' -gCg 3'	5' -gCg 3'	5' -gCg 3'	5' -gCg 3'
3'-primer³	87(346)	86(344)	86(344)	40(206)	70(296)
	5' -CTC 3'	5' -CAC 3'	5' -CCA 3'	5' -CgA 3'	5' -TCT 3'
Well No.	1	2	3	4	5
DRB1 allele⁴					
*10:01:01-10:01:03	1	2			
*10:02	1		3		
*10:03	1	2		4	
*10:04	1	2			5
DRB1 allele					
Well No.	1	2	3	4	5