

Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*12:05	1	2			5	6	7	8				
*12:06	1	2		4	5	6	7	8	9			
*12:07	1	2		4	5	6	7		9			
*12:08	1	2		4	5	6		8	9			
*12:09	1	2	3	4		6	7	8	9			
*12:10	1	2		4	5	6	7	8	9			
*12:11	1	2		4	5	6	7	8	9			
*12:12	1	2		4	5	6	7		9			
*12:13	1	2		4	5	6	w		9	10		
*12:14	1	2			5	6	7	8				
*12:15	1	2			5	6	7			10		
*12:16:01	1	2		4	5					10		
*12:16:02	1	2		4	5				9	10		
*12:17	1	2		4	5	6	7	8	9			
*12:18	1	2		4		6	7		9	10		
*12:19	1	2		4	5				9	10	11	
*12:20	1	2		4	5	6	7		9	10		
*12:21		2		4	5	6	7		9	10		
*12:22	1			4	5			8				
*12:23	1	2		4	5	6			9	10		
*12:24N	1	2		4	5	6	7	8	9			
*12:25	1	2		4	5	6	7	8	9			
*12:26	1	2		4	5	6	7		9	10		
*12:27	1	2		4	5	6			9	10		
*12:29, 12:36	1	2		4	5	6	7	8	9			
*12:31N	1	2		4	5	6	7		9	10		
*12:33	1	2		4		6	7		9	10		
*12:37	1	2		4	5	6			9	10		
*12:38		2			5	6	7	8				
*12:39	1	2		4	5			8	9			
*01:01:01-01:59												
*03:92	1		3								11	
*08:01:01-08:01:05, 08:07-08:08, 08:11, 08:16, 08:26, 08:39, 08:43, 08:50										10		
*08:02:01-08:02:04, 08:09, 08:21, 08:42, 08:44			3							10		
*08:03:02-08:03:03, 08:14-08:15, 08:23, 08:27, 08:29, 08:30:02, 08:33, 08:35-08:36:02, 08:38, 08:46, 08:49, 08:51								8				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
Length of spec.	135	215	165	105	165	250	215	195	165	195	250	170
PCR product(s)												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12
*08:04:01, 08:04:04-08:04:05, 08:04:07, 14:15			3							10	11	

13	14	15	16	17	18	19	20	21	22	23	24	Well No.
13								21				*12:05
	14							21				*12:06
		15						21				*12:07
			16					21				*12:08
								21				*12:09
				17				21				*12:10
					18			21				*12:11
						19		21				*12:12
							20	21				*12:13
13												*12:14
13								21				*12:15
								21	22			*12:16:01
								21	22			*12:16:02
								21		23		*12:17
								21			24	*12:18
								21				*12:19
13								21				*12:20
								21	22			*12:21
								21	22			*12:22
			16				20	21				*12:23
13								21				*12:24N
				17				21		23		*12:25
					18			21				*12:26
			16					21				*12:27
												*12:29, 12:36
			16					21				*12:31N
								21				*12:33
								21				*12:37
									22			*12:38
								21	22			*12:39
										23		*01:01:01-01:59
												*03:92
												*08:01:01-08:01:05, 08:07-08:08, 08:11, 08:16, 08:26, 08:39, 08:43, 08:50
												*08:02:01-08:02:04, 08:09, 08:21, 08:42, 08:44
												*08:03:02-08:03:03, 08:14-08:15, 08:23, 08:27, 08:29, 08:30:02, 08:33, 08:35-08:36:02, 08:38, 08:46, 08:49, 08:51
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
115	135	200	75	90	90	195	220	115	105	110	170	Length of spec.
185			115		135				220			PCR product(s)
255												
13	14	15	16	17	18	19	20	21	22	23	24	Well No.
												*08:04:01, 08:04:04-08:04:05, 08:04:07, 14:15

