

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
<b>DRB1 allele<sup>1,2,3</sup></b>																								
*12:01:01, 12:01:04-12:01:05, 12:28, 12:30, 12:34-12:35, 12:40-12:41	1	2		4	5	6	7	8	9												21			
*12:01:02-12:01:03	1	2		4	5	6	7	8													21			
*12:02:01, 12:02:03-12:02:06, 12:32, 12:44	1	2		4	5	6	7		9	10											21			
*12:02:02	1	2		4	5	6	7			10											21			
*12:03:02	1	2		4	5			8			11										21			
*12:03:03	1	2		4	5			8	9		11										21			
*12:04	1	2		4		6	7	8	9			12									21			
*12:05	1	2			5	6	7	8					13								21			
*12:06	1	2		4	5	6	7	8	9					14							21			
*12:07	1	2		4	5	6	7		9						15						21			
*12:08	1	2		4	5	6		8	9							16					21			
*12:09	1	2	3	4		6	7	8	9												21			
*12:10	1	2		4	5	6	7	8	9								17				21			
*12:11	1	2		4	5	6	7	8	9									18			21			
*12:12	1	2		4	5	6	7		9										19		21			
*12:13	1	2		4	5	6	w		9	10										20	21			
*12:14	1	2			5	6	7	8					13											
*12:15	1	2			5	6	7			10			13								21			
*12:16:01, 12:16:03	1	2		4	5					10											21	22		
*12:16:02	1	2		4	5				9	10											21	22		
*12:17	1	2		4	5	6	7	8	9												21		23	
*12:18	1	2		4		6	7		9	10					15						21			
*12:19	1	2		4	5				9	10	11										21			
*12:20	1	2		4	5	6	7		9	10			13								21			
*12:21		2		4	5	6	7		9	10											21	22		
*12:22	1			4	5			8													21	22		
*12:23	1	2		4	5	6			9	10						16				20	21			
*12:24N	1	2		4	5	6	7	8	9				13								21			
*12:25	1	2		4	5	6	7	8	9								17				21		23	
*12:26	1	2		4	5	6	7		9	10								18			21			
*12:27	1	2		4	5	6			9	10						16					21			
*12:29, 12:36	1	2		4	5	6	7	8	9															
*12:31N	1	2		4	5	6	7		9	10						16					21			
*12:33	1	2		4		6	7		9	10											21			
*12:37	1	2		4	5	6			9	10											21			
*12:38		2			5	6	7	8														22		
*12:39	1	2		4	5			8	9												21	22		
*12:42							7		9															
*12:43	1			4	5	6	7			10											21			
*12:45	1	2		4	5	6	7		9	10														
*12:46	1	2		4		6	7	8	9												21			
*12:47	1	2		4		6	7	8	9						15						21			
*12:48	1	2	w	4		6	7	8	9												21			
*12:49	1	2		4		6	7		9	10		12									21			
*01:01:01-01:63																							23	
*03:92	1		3								11													
*08:01:01-08:01:05, 08:07-08:08, 08:11, 08:16, 08:26, 08:39, 08:43, 08:50, 08:55										10														
*08:02:01-08:02:04, 08:09, 08:21, 08:42, 08:44				3						10														
*08:03:02-08:03:04, 08:15, 08:23, 08:27, 08:29, 08:30:02, 08:33, 08:35-08:36:02, 08:38, 08:46, 08:49, 08:51, 08:56, 08:58								8																

Negative Control

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>
*08:04:01, 08:04:04-08:04:05, 08:04:07, 14:15			3							10	11													
*08:04:02-08:04:03			3							10	w													
*08:04:06, 08:06, 08:54, 14:73										10	11													
*08:05		2								10														
*08:10								8			11													
*08:12					6			8																
*08:13			3																19					
*08:17	1									10														
*08:18, 08:40, 08:47		2						8																
*08:19, 08:34					5			8																
*08:22						6				10														
*08:24		2	3							10														
*08:25		2			5			8																
*08:28	1		3							10	11													
*08:30:01, 08:30:03			3					8																
*08:31		2								10	11	12												
*08:32				4				8														21	22	
*08:37	1							8																
*08:41		2								10		12												
*08:45	1		3					8																
*08:48																			19					
*08:52					5					10														
*08:53				4				8								16						21		
*11:67	1									10	11	12												
*11:76, 11:114, 13:34, 13:64, 13:136, 13:174, 14:41, 14:77, 14:110																16								
*13:17, 13:175	1	2	3								11													
*13:116, 14:52		2	3								11													
*13:145																							22	
*14:04, 14:50, 14:76, 14:79, 14:107, 14:120, 14:145, 14:152N											11													
*14:11, 14:148											11	12												
*14:28						6																		
*14:31		2									11													
*14:126			3								11													
*14:138	1										11													
*15:01:01:01-15:01:16, 15:01:18- 15:101, 15:103-15:112, 16:01:01- 16:05:02, 16:07-16:24															w									
*15:102															w	16								
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>

Negative Control