

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DRB1 allele^{1,2,3}																								
*12:01:01:01-12:01:01:06, 12:01:04-12:01:09, 12:28, 12:30, 12:34-12:35, 12:40-12:41, 12:54, 12:59, 12:68, 12:70, 12:72N, 12:77, 12:79	1	2		4	5	6	7	8	9												21			
*12:01:02-12:01:03	1	2		4	5	6	7	8													21			
*12:02:01:01-12:02:01:04, 12:02:03-12:02:06, 12:02:08-12:02:09, 12:32, 12:44, 12:50, 12:61, 12:67, 12:69, 12:71, 12:74N-12:75, 12:78, 12:80	1	2		4	5	6	7		9	10											21			
*12:02:02, 12:55	1	2		4	5	6	7			10											21			
*12:02:07, 12:51-12:52	1	2		4	5	6	7		9	10														
*12:03:02	1	2		4	5			8			11										21			
*12:03:03	1	2		4	5			8	9		11										21			
*12:04	1	2		4		6	7	8	9			12									21			
*12:05	1	2			5	6	7	8					13								21			
*12:06	1	2		4	5	6	7	8	9					14							21			
*12:07	1	2		4	5	6	7		9						15						21			
*12:08	1	2		4	5	6		8	9							16					21			
*12:09	1	2	3	4		6	7	8	9												21			
*12:10	1	2		4	5	6	7	8	9								17				21			
*12:11	1	2		4	5	6	7	8	9									18			21			
*12:12	1	2		4	5	6	7		9										19		21			
*12:13	1	2		4	5	6	w		9	10										20	21			
*12:14	1	2			5	6	7	8					13											
*12:15	1	2			5	6	7			10			13								21			
*12:16:01, 12:16:03	1	2		4	5					10											21	22		
*12:16:02	1	2		4	5				9	10											21	22		
*12:17	1	2		4	5	6	7	8	9												21		23	
*12:18	1	2		4		6	7		9	10					15						21			
*12:19	1	2		4	5				9	10	11										21			
*12:20	1	2		4	5	6	7		9	10			13								21			
*12:21		2		4	5	6	7		9	10											21	22		
*12:22	1			4	5			8													21	22		
*12:23	1	2		4	5	6			9	10						16				20	21			
*12:24N	1	2		4	5	6	7	8	9				13								21			
*12:25	1	2		4	5	6	7	8	9				13				17				21			
*12:26	1	2		4	5	6	7		9	10								18			21			
*12:27	1	2		4	5	6			9	10						16					21			
*12:29, 12:36, 12:63	1	2		4	5	6	7	8	9															
*12:31N	1	2		4	5	6	7		9	10						16					21			
*12:33	1	2		4		6	7		9	10											21			
*12:37, 12:64-12:65	1	2		4	5	6			9	10											21			
*12:38		2			5	6	7	8														22		
*12:39	1	2		4	5			8	9												21	22		
*12:42							7		9															
*12:43	1			4	5	6	7			10											21			
*12:45	1	2		4	5	6	7		9	10									19					
*12:46	1	2		4		6	7	8	9												21			
*12:47	1	2		4		6	7	8	9					15							21			
*12:48	1	2	3	4		6	7	8	9								17				21			
*12:49	1	2		4		6	7		9	10		12									21			
*12:53		2		4	5	6	7	8	9												21	22		
*12:56	1	2			5	6	7			10			13											
*12:57																				20				
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24

Negative Control

