

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
DRB1 allele <sup>1,2,3</sup>																
*07:01:01:01-07:01:08, 07:29-07:33, 07:35	1	2														
*07:03	1		3													
*07:04	1			4												
*07:05	1	2			5											
*07:06	1	2				6										
*07:07	1	2					7									
*07:08	1	2						8								
*07:09	1	2	3													
*07:10N	1	2		4												
*07:11	1				5											
*07:12	1	2							9							
*07:13	1	2								10						
*07:14, 07:24 <sup>4</sup>	1	2									11					
*07:15	1	2										12				
*07:16	1	2											13			
*07:17, 07:27 <sup>4</sup>	1	2												14		
*07:18	1	2								10						15
*07:19	1	2														15
*07:20	1	2												14	15	
*07:21	1	2										12	13			
*07:22		2							9							
*07:23	1	2							9			12				
*07:25				4	5											
*07:26N	1	2	3					8								
*07:28	1	2					7									15
*07:34	1	2								10	11					
*04:90, 09:08									9							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Negative Control