

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
DRB1 allele^{1,2,3}																								
*08:01:01-08:01:02, 08:01:04-08:01:06, 08:50, 08:77, 08:79	1		3		5		7																	
*08:02:01:01-08:02:04, 08:42, 08:72-08:73, 08:82		2	3		5		7																	
*08:03:02:01-08:03:07, 08:46, 08:51, 08:53, 08:58, 08:61-08:62, 08:66, 08:68:01-08:69, 08:71, 08:74, 08:76	1			4	5		7									16								
*08:04:01, 08:04:04-08:04:07, 08:59		2	3		5		8										17					23		
*08:04:02-08:04:03		2	3		5		8										w							
*08:05	1		3			6	7																	
*08:06, 08:54, 08:70	1		3		5		8										17			20			23	
*08:07			3		5		7												19					
*08:08			3		5		7			11														
*08:09		2	3		5		7											18						
*08:10	1			4	5		8									16	17			20			23	
*08:11			3		5		7						13											
*08:12	1			4	5		8		10							16	17			20				
*08:13		2			5		7				12													
*08:14	1			4			7		9							16	17							
*08:15				4	5		7				11													
*08:16	1		3		5		7													19				
*08:17	1		3		5		7							14								21		
*08:18, 08:47	1			4		6	7									16								
*08:19				4	5		7								15								22	
*08:20, 14:12:01-14:12:02, 14:84, 14:156, 14:173 ⁴					5												17	18						
*08:21		2	3				7											18						
*08:22	1		3		5		8		10								17			20				
*08:23	1			4	5		7									16							22	
*08:24:01-08:24:02		2	3			6	7																	
*08:25				4		6	7								15									
*08:26	1		3		5		7		9															
*08:27	1			4	5		7					12				16								
*08:28		2	3		5		8						14				17					21	23	
*08:29	1			4			7									16						22		
*08:30:01-08:30:03		2		4	5		7																	
*08:31			3			6	8			11													23	
*08:32	1			4	5		7									16		18						
*08:33	1			4	5		7						13			16								
*08:34				4	5		7								15									
*08:35	1			4	5		7		9							16		18						
*08:36:01-08:36:02	1			4	5		7		9							16								
*08:37	1			4	5		7							14		16						21		
*08:38	1			4	5		7									16			19					
*08:39	1		3		5		7						13											
*08:40	1			4		6	7		10							16								
*08:41			3			6	7				11													
*08:43	1		3				7					12											23	
*08:44		2	3		5		7					12												
*08:45:01-08:45:02		2		4	5		7							14								21		
*08:48	1				5		7					12												
*08:49	1			4	5		7								15	16			19					
*08:52			3		5		7								15									
*08:55, 08:64	1		3		5																			
*08:56	1			4			7									16								
*08:57					5																			
*08:60N	1			4	5		7							14		16								
*08:63				4	5		7				11			14								21		
*08:65	1			4	5											16								
*08:67		2	3				8										17						23	
*08:75			3		5		8			11							17						23	
*08:78N	1		3		5		7															21		
*08:80					5												17							
*08:81		2			5		7																	
*03:12, 14:65, 14:169-14:170									9													20		
*03:92		2					8							14									23	
*04:10:01-04:11:04, 04:67, 04:91, 04:126, 04:147, 04:201, 04:203, 04:207, 13:75, 13:108, 13:169, 15:12																						20		
*04:12																16	17				20			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*04:18, 04:25, 04:58																	17							

Negative Control

