

*04:10:01-04:11:04, 04:11:06, 04:67, 04:91, 04:126, 04:147, 04:201, 04:203, 04:207, 04:250, 13:75, 13:108, 13:169, 13:276, 15:12																			20					
*04:12														16	17			20						
*04:18, 04:25, 04:58, 11:261															17									
*04:86, 04:106, 13:03:01-13:03:09, 13:12:01-13:12:04, 13:30, 13:33:01-13:33:03, 13:38, 13:66:01-13:66:02, 13:95, 13:101, 13:115, 13:118, 13:122, 13:151-13:152, 13:161, 13:167, 13:174, 13:188, 13:194, 13:216, 13:219, 13:227-13:228, 13:253, 13:295N, DRB3*02:131														16										
*04:136, 13:04, 13:48, 13:58, 13:81, 13:89:01-13:89:02, 13:94:01-13:94:02														16					20					
*04:202, 14:13									9															
*11:25, 11:119, 13:18, 13:144, 13:156					5												17	18			21			
*11:67			3		5			8		11			14			17					21		23	
*11:193:01-11:193:02			3			6	7			11			14											
*11:209			3			6	8			11			14										23	
*11:229						6	8		10	11			14										23	
*12:01:01:01-12:01:07, 12:01:09, 12:05-12:06, 12:08, 12:10-12:11, 12:14, 12:17, 12:24N-12:25, 12:28-12:30, 12:34-12:36, 12:40-12:41, 12:54, 12:58-12:59, 12:63, 12:68, 12:70, 12:72N, 12:77, 12:79, 12:82				4		6	8		10				14	15										
*12:01:08, 12:47				4		6	8		10				14											
*12:02:01:01-12:02:05, 12:02:07-12:02:10, 12:13, 12:15, 12:20, 12:23, 12:26-12:27, 12:31N-12:32, 12:37, 12:44-12:45, 12:50-12:52, 12:55-12:56, 12:60N-12:61, 12:64-12:65, 12:67, 12:69, 12:71, 12:74N-12:75, 12:78, 12:80-12:81				3		6	8		10				14	15										
*12:02:06, 12:18, 12:33				3		6	8		10				14											
*12:03:02-12:03:03				4		6	8		10				14	15									23	
*12:04				4		6	8		10	11			14											
*12:07, 12:73						6	8		10				14	15										
*12:09			2		4	6	8		10				14											
*12:12						6	8		10		12		14	15										
*12:16:01-12:16:03				3		6	7						14	15										
*12:19				3		6	8						14	15									23	
*12:21				3		6	8		10				15											
*12:22					4		7						14	15										
*12:38, 12:53				4		6	8		10				15											
*12:39				4		6	7						14	15										
*12:43				3			8		10				14	15										
*12:46		1		4		6	8		10				14	16							20			
*12:48			w	4		6	8		10				14											
*12:49				3		6	8		10	11			14											
*12:62			2			6	8		10		12		14											
*12:66							8		10				14	15									22	
*12:76			2	3		6	8		10															
*13:13					5									16		18								
*13:17, 13:175, 13:214			2		4	6	8		10				14										23	
*13:32, 13:93									9					16									20	
*13:47, 13:55, 13:232, 14:03:01-14:03:02, 14:27:01-14:27:02, 14:40, 14:55, 14:67, 14:77, 14:89, 14:102, 14:115-14:116, 14:135-14:136, 14:144, 14:164, 14:174, 14:181, 14:218					5										18									
*13:65, 13:120, 13:139									9					16										
*13:116			2		4	6	8		10														23	
*13:154					5									16		18						21		
*13:164					5										17	18						20		
*14:04:01:01-14:04:07, 14:11, 14:120, 14:145, 14:148, 14:152N, 14:180, 14:196, 14:201, 14:205-14:206, 14:210Q							8		11														22	23
*14:15			2	3		5		8							17	18								23
*14:28							8		10	11													22	
*14:31						6	8			11													23	
*14:50, 14:76, 14:79, 14:107							8																23	
*14:52			2			6	8																23	
*14:63						5			9					16		18								
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
*14:68:01-14:68:02, 14:93							7			11													22	
*14:71								8		11													22	
*14:73				3				8		11													22	23
*14:78					5				9						16	17	18				20			
*14:85					5				9								18							

Negative Control

e Control

*14:126:01-14:126:02		2						8													22	23	Negati	
*14:138								8			11			14							22	23		
*14:175		2						8						14							22	23		
*15:21																17					21			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24