

Well No.	1	2	3	4	5	6	7	8	9	10
DRB1 allele^{1,2,3}										
*09:01:02:01-09:01:14, 09:26, 09:28-09:30, 09:33-09:41, 09:43, 09:45-09:52N	1	2								
*09:02:01-09:02:02	1		3							
*09:03, 09:15	1	2		4						
*09:04	1	2			5					
*09:05		2				6				
*09:06	1	2					7			
*09:07		2						8		
*09:08	1				5					
*09:09, 09:13, 09:17-09:18 ⁴	1	2								
*09:10-09:11, 09:14	1	2				6				
*09:12, 09:19	1	2						8		
*09:16		2		4						
*09:20		2						8		
*09:21:01:01-09:21:01:02	1	2								10
*09:22	1	2							9	
*09:23	1									
*09:24, 09:32, 09:44	1									
*09:25	1	2								
*09:27	1	2								
*09:31:01-09:31:02	1	2								
*09:42	1		3		5					
*09:53		2		4		6				
*01:65:01, DRB4*01:125					5	6				
*01:136, DRB4*01:124										10
*04:90, 04:155, 04:167, 07:50, 08:108, 13:220, 16:42, DRB5*01:12, DRB5*01:15, DRB5*02:31, DRB5*02:34										
*07:52								8		

11	12	13	14	15	16
11					
		14			
11					
	12				
		13			
				15	
11					