

INTERPRETATION TABLE**DRB1*1401/1454 SSP typing**

	Well	
	1	2
Length of spec.	250	210
PCR product		
Length of int.	515	430
pos. control ¹		
5'-primer(s) ²	112(421)	13(125)
	5' -ACT 3'	5' -gTC 3'
3'-primer(s) ³	181(630)	70(296)
	5' -CTT 3'	5' -TCC 3'
Well No.	1	2
DRB1 allele		
*140101-140103	1	2
*1454		2
*140501-140503, 140701-1408, 1414, 1418, 142301-142303, 1426, 143201, 143202, 1434-1436, 1438, 1439, 1442-1445, 1456, 1458-1460, 1462, 1464, 1465, 1470, 1472, 1475, 1481, 1482, 1486-1488, 1490-1492N, 111301, 111302, 1117, 1152		2
*1455		w
DRB1 allele		
Well No.	1	2