

| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 |
|---|----------|----------|----------|----------|----------|----------|----------|----------|----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| DRB5 allele^{1,2,3} | | | | | | | | | | | | | | | | | | | | | |
| *01:01:01:01-01:01:01:04, 01:01:04-01:01:06, | 1 | 2 | 3 | 4 | 5 | | | | | | | | | | | | | | | | |
| *01:01:02 | 1 | 2 | 3 | 4 | | | | | | | | | ? | | | | | | | | |
| *01:01:03, 01:66, 01:75, 01:86 | 1 | 2 | 3 | 4 | | | | | | | | | | | | | | | | | |
| *01:02, 01:25, 01:28, 01:32, 01:35, 01:59- 01:60, 01:64, 01:67N-01:68N, 01:70, 01:72, 01:74, 01:76 | 1 | 2 | 3 | | | 6 | 7 | | | | | | | | | | | | | | |
| *01:03 | 1 | 2 | | | | 6 | 7 | 8 | | | | | ? | | | | | | | | |
| *01:04 | 1 | 2 | 3 | 4 | | | | | 9 | | | | | | | | | | | | |
| *01:05 | 1 | 2 | 3 | | 5 | 6 | | | | | | | | | | | | | | | |
| *01:06 | | | | 4 | | | | 8 | | | 11 | | | | | | | | | | |
| *01:07 | 1 | 2 | 3 | 4 | 5 | | | | | | | | | | | | | | | | |
| *01:08:01N-01:08:02N, 01:56 | 1 | 2 | 3 | | | 6 | 7 | | | | | | 13 | | | | | | | | |
| *01:09 | 1 | 2 | 3 | 4 | 5 | | | 8 | | | | | ? | | | | | | | | |
| *01:10N | 1 | 2 | 3 | | | 6 | 7 | | | | | | | | 15 | | | | | | |
| *01:11 | 1 | | | 4 | | | | 8 | | | 11 | | | | | | | | | | |
| *01:12 | 1 | 2 | 3 | | | | | | | 10 | | | | | 15 | | | | | | |
| *01:13, 01:41 | 1 | 2 | 3 | | 5 | | | | 9 | | | | | | | | | | | | |
| *01:14 | 1 | 2 | 3 | | | | | | | | 11 | | | | | | | | | | |
| *01:15 | 1 | 2 | 3 | 4 | | | | | | 10 | | | | | 15 | | | | | | |
| *01:16 | 1 | 2 | 3 | 4 | 5 | | | | 9 | | | | | | | | | | | | |
| *01:17 | 1 | 2 | 3 | 4 | 5 | | 7 | | | | | | | | | | | | | | |
| *01:18 | 1 | 2 | 3 | 4 | 5 | | | | | | | | 13 | | | | | | | | |
| *01:20 | 1 | 2 | 3 | | | 6 | 7 | | | | | | | | 16 | | | | | | |
| *01:21 | | | | 4 | | | | 8 | | | 11 | | | | 16 | | | | | | |
| *01:27N, 01:46 | 1 | 2 | 3 | | | 6 | 7 | | | | | | | | | | 17 | | | | |
| *01:31, 01:44 | | 2 | 3 | 4 | 5 | | | | | | | | | | | | | | | | |
| *01:39 | | | | | | 6 | 7 | | | | | | | | | | | | | | |
| *01:43 | | 2 | | 4 | | | | 8 | | | | | | | | 16 | | | | | |
| *01:47 | 1 | 2 | 3 | 4 | 5 | | | | | | | | | | 16 | | | | | | |
| *01:48N | 1 | 2 | 3 | 4 | 5 | | | | | 10 | | | | | | | | | | | 20 |
| *01:49N | 1 | 2 | 3 | 4 | 5 | | | | | | | | | | | | | | | 19 | |
| *01:51, 01:63:02, 01:73, 01:78, 01:87 | 1 | 2 | 3 | | 5 | | | | | | | | | | | | | | | | |
| *01:52N, 01:77 | 1 | 2 | 3 | | | | | | | | | | | | | | | | | | |
| *01:55 | | 2 | | 4 | | | | 8 | | | | | | | | | | | | | |
| *01:61 | 1 | 2 | 3 | | | | 7 | | | | | | | | | | | | | | |
| *01:69 | 1 | 2 | | | | 6 | 7 | 8 | | | | | | | | | 17 | | | | |
| *01:71N | 1 | | | 4 | | | | | | | | | | | | | | | | | |
| *01:79Q, DRB1*16:60 | | | | 4 | 5 | | | | | | | | | | | | | | | | |
| *01:82 | 1 | 2 | | 4 | 5 | | | | | | | | | | | | | | | | |
| *02:02:01-02:02:03, 02:09, 02:11, 02:14- 02:16, 02:18, 02:20-02:23, 02:27-02:29 | | | | | | | | 8 | | | 11 | 12 | | | | | | | | | |
| *02:03 | 1 | | | | | | | 8 | | | 11 | | | | | | | | | | |
| *02:04 | | 2 | | | | | | 8 | | | | 12 | ? | | | | | | | | |
| *02:05 | | | 3 | | | | 7 | | | | | 12 | | | | | | | | | |
| *02:06 | | | | 4 | | | | 8 | | | 11 | 12 | | | | | | | | | |
| *02:07 | | | | | | | | 8 | | | 11 | 12 | | 14 | | | | | | | |
| *02:08 | | 2 | 3 | | | 6 | 7 | | | | | 12 | | | 16 | | | | | | |
| *02:10 | | | | | | | | 8 | | | 11 | 12 | | | 16 | | | | | | |
| *02:12 | | 2 | 3 | | | 6 | 7 | | | | | 12 | | | | | 17 | | | | |
| *02:13 | | | | | | | | 8 | | | 11 | | | | | | | | | 18 | |
| *02:17 | | | | | | | | 8 | | | 11 | | | | | | | | | | |
| *02:19N | | | | | | | | 8 | | | 11 | 12 | | | | | 17 | 18 | | | |
| *02:24 | | | | | | | | 8 | | | | 12 | | | | | | | | | |
| *02:25N | | 2 | 3 | | | 6 | 7 | | | | | 12 | | | | | | | | | |
| *02:26N | | 2 | 3 | | | 6 | 7 | | | | | 12 | 13 | | | | | | | | |
| DRB1*03:112, DRB1*09:06, DRB4*01:117 DRB1*08:59, DRB1*11:210 | | | | | | | | | | | | | | | | 16 | | | | 18 | |
| DRB1*09:07 | 1 | 2 | 3 | | | | | | | | | | | | 15 | | | | | | |
| DRB1*15:02:03, DRB1*16:01:02, DRB1*16:02:02, DRB1*16:05:01, DRB4*01:05, DRB4*01:07:01 DRB1*15:86 | | | | | | | | | | | | | ? | | | | | | | | |
| | | | | | | | | | | | | | 13 | | | | | | | | |
| Well No. | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 |

Negative Control