

Tube No.		1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
HLA-A allele ¹	ser ²																				
*01010101-0104N, 0106, 0107, 0109-0111N, 0116N-0118N, 0120, 0122N-0127N	A1, null	1			4																
*0108, 0114, 0115N	1	1																			
*0112, 0119	-	1		3																	
*0113	-	1			4							11		13							
*0121	-	1		3	4																
*02010101-0209, 0211-0216, 0218-0222, 0224-0226, 0228-0233, 0236-0238, 0240-0245, 0247, 0249, 0251, 0253N, 0254, 0257-0261, 0263-0269, 0271, 0272, 027401-0277, 0279-0292, 0294N, 0296, 0297, 0299, 9201, 9202, 9204-9209, 9211-9213N, 9215, 9216, 9218-9221, 9	A2, A203, null		2																		
*0210, 0239, 0250, 0252, 0273, 0293, 0295, 9210, 9214, 9217	A210		2																		20
*021701-021702	A2		2					w													
*0227	-		2									11									
*0234, 023501, 023503, 0256, 0262, 9203	-		2											13							20
*023502	-													13							20
*0246, 0270	-		2																		20
*0248	-																				20
*0255	-		2						8												
*0278	-		2		4																20
*9222	-						6														20
*9224	-		2							9											
*03010101-0307, 0309-0311N, 0313-0317, 0319-0323, 0325-0329	A3			3										13							
*0308	-			3																	
*0312	-			3	4																
*0318, 2619 ³	-				4									13							
*0324	-			3					8					13							
*0330	-			3		5								13							
*110101-1109, 1112-1124, 1126-1132	A11				4							11									
*1110	-				4				8			11									
*1111	-				4							11	12								
*1125	-			3	4							11									
*2301-2313, 2315, 2316	A23					5	6														
*2314, 2405, 2424, 2465 ⁴	-					5	6	7													
*24020101-2404, 2406-2411N, 2413, 2417, 2418, 2420-2423, 2425-2443, 2445N-2450, 2454-2456, 2458-2463, A24 2466-2479						5		7													
*2414, 2415, 2451-2453, 2457, 2464	A24					5															
*2419, 2444	-					5		7				11									
*250101-2504	A25								8	9				13							
*2505	-								8	9				13							20
*2506	-								8	9			12								
*260101-260105, 2608, 2610-2615, 2617, 2618, 2620, 2623-2629, 2632, 2633, 2635	A26, null								8		10			13							
*2602	A26								8		10			w							
*2603, 2606, 2621	A26								8			11		13							
*2604	A26								8		10										
*2605, 2630	A26								8					13							
*260701-260702	A26										10			13							
*2609	A26								8		10		12								
*2616	-							7			10			13							
*2622	-								8		10			13							
*2631	-								8					13							
*2634	-								8		10										
*29010101-2906, 2908N-2912, 2915, 2916	A29														14						
*2907	-						6								14						

