

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	
*66:02, 66:16, 66:21, 66:25-66:26Q, 66:34, 66:43								8				12												24				28					
*66:03:01:01-66:03:01:03, 66:28N, 66:39N								8				12											24										
*66:06								8			11	13									21										30		
*66:09								8			11	13	14								21										30		
*66:10								8			11	w																		30			
*66:12								8				13											24				28			30			
*66:19								8			11	12	13																		30		
*66:23				4							11	13																			30		
*66:42											11	13											24								30		
*68:01:01:01-68:01:17, 68:01:19-68:02:16, 68:06-68:08:01:02, 68:09-68:10, 68:12:01-68:12:02, 68:13:02-68:14, 68:16-68:18N, 68:21:01-68:24, 68:26-68:28, 68:30, 68:32-68:40, 68:42-68:44, 68:46-68:65, 68:67-68:70, 68:73-68:83, 68:85-68:89, 68:91-68:92, 68:94N-68:96, 68:98-68:102, 68:104:01-68:104:02, 68:106-68:108, 68:110-68:116, 68:118-68:129, 68:131-68:133, 68:135, 68:137-68:156, 68:158-68:177, 68:179, 68:181N-68:205N, 68:207-68:211, 68:213N-68:217, 68:219-68:226, 68:228-68:229, 68:231-68:240	A28, A68, Null, -							8													20								28				
*68:01:18, 68:03:01:01-68:04:02, 68:31, 68:72, 68:90, 68:93, 68:97, 68:180, 68:206, 68:212	A68, -							8													20												
*68:05, 68:15, 68:20	A68							8													20		23							29			
*68:08:02, 68:105, 68:157	A68, -							8													20	21							28				
*68:11N, 68:41, 68:230	Null, -							8													20	22							28				
*68:13:01, 68:66:01-68:66:02, 68:134	A68, -			4				8													20								28				
*68:19, 68:130:02	A68, -							8				14									20								28				
*68:25	A68							8												19	20							28					
*68:29	A68							8												18	20	21						28					
*68:45, 68:117	-					7		8													20								28				
*68:71	-							8				13									20								28				
*68:84	-							8		10											20								28				
*68:103:01-68:103:02	-							8													20							27	28				
*68:109	-							8													20		23										
*68:130:01	-							8				12									20					26		28					
*68:136	-							8													20		23					28					
*68:178	-			4																	20								28				
*68:218	-							8																					28				
*68:227	-							8		10											20												
*69:01:01:01-69:01:03, 69:03-69:06	A69, -							8																21				28					
*69:02	-							8			11												21										
*69:07	-							8					13										21					28					
*69:08	-							8															21										
*74:01:01:01-74:01:02, 74:01:04-74:03:02, 74:05, 74:08-74:12N, 74:14N-74:18, 74:20, 74:22, 74:24-74:29, 74:31-74:32N, 74:34, 74:36, 74:38	A74, Null, -																				19				25								
*74:01:03	-											12											19	24	25								
*74:04	A74																				18	19			25								
*74:06, 74:39	A74, -																					19			25				29				
*74:07, 74:33	A74, -																				17	19			25								
*74:13	-												13	16								19			25		28						
*74:19	-			4																		19			25								
*74:21	-																				18	19			25				29				
*74:23	-																					19			25		27						
*74:30	-											11										19			25								
*74:35	-																					19		21		25							
*74:37	-												13									19			25								

Negative Control

