

INTERPRETATION TABLE

HLA-A*11 Add-on SSP typing

	Well		
	1	2	3
Length of spec.	270	210	170
PCR product			
Length of int.	800	1070	1070
pos. control ¹			
5'-primer ²	28	831	874
	5' -TCC 3'	5' -gAg 3'	5' -CCg 3'
3'-primer ³	127	899	899
	5' -CTT 3'	5' -ACg 3'	5' -ACA 3'
Well No.	1	2	3
HLA-A allele			
*11:02:01-11:02:03, 11:14, 11:16, 11:38, 11:57, 11:101	1		
*11:77	1		3
*11:110	1	2	
*02:294, 66:08		2	
*03:26			3
HLA-A allele			
Well No.	1	2	3