

INTERPRETATION TABLE		
HLA-A*11 Add-on SSP typing		
	Well	
	1	2
Length of spec.	270	210
PCR product		
Length of int.	800	1070
pos. control¹		
5'-primer²	28	831
	5' -TCC 3'	5' -gAg 3'
3'-primer³	127	899
	5' -CTT 3'	5' -ACg 3'
Well No.	1	2
HLA-A allele		
*11:02:01-11:02:03, 11:14, 11:16, 11:38, 11:57, 11:77, 11:101, 11:113	1	
*11:110	1	2
*02:294, 32:54, 66:08		2
HLA-A allele		
Well No.	1	2