



























<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>*02:01:01:01</b> , 02:01:01:03- <b>02:01:02</b> , <b>02:01:04</b> -02:01:15, <b>02:01:18</b> -02:01:19, 02:01:21- <b>02:01:51</b> , 02:01:53-02:01:62, 02:01:64-02:01:73, 02:01:75- 02:01:81, <b>02:01:84</b> -02:01:86, 02:01:89-02:01:90, 02:01:92- 02:01:99, 02:201, 02:204, 02:208, 02:210, 02:212, 02:216-02:217:01, 02:218, 02:220, 02:231, 02:234, 02:236, 02:238-02:241, 02:249, 02:252, 02:256- 02:257, 02:262, 02:266, 02:270, 02:272-02:273, 02:276 02:277, 02:285, 02:287- 02:288, 02:292, 02:294, 02:296, 02:302, 02:306- 02:307, 02:311-02:313, 02:316, 02:318, 02:326- 02:327, 02:329, 02:332, 02:336, 02:340-02:341, 02:343, 02:346, 02:348- 02:349, 02:352-02:354, 02:357, 02:360, 02:362- 02:363, 02:365, 02:368, 02:374-02:375, 02:379-02:381, 02:383, <b>02:384</b> , 02:385- 02:386, 02:388-02:392, 02:396 02:397, 02:401, 02:410- 02:411, 02:416, 02:423- 02:424, 02:430, 02:434- 02:435, 02:441-02:446, 02:448, 02:456 <sup>5</sup>
<b>*02:01:01:02L</b>
<b>*02:01:03</b> , 02:243:01- 02:243:02, 02:394
*02:01:17, <b>02:01:52</b> , 02:01:63, 02:01:74, 02:372, 02:378, 02:414
*02:01:83, 02:399
*02:01:87, 02:01:91, 02:317, 02:347, 02:361, 02:422
*02:01:88, 02:418
*02:02, 02:209, 02:320
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>

*02:03:01, 02:03:03-02:03:06, 02:253, 02:258, 02:264, 02:281, 02:345, 02:370, 02:427, 02:431
*02:03:02
*02:04
*02:05:01-02:05:06, 02:229, 02:232, 02:286, 02:337, 02:359
*02:06:01-02:06:11, 02:06:13- 02:06:15, 02:248, 02:278, 02:290, 02:328, 02:333, 02:404-02:405, 02:428
*02:06:12
*02:07:01-02:07:03, 02:219, 02:261, 02:265, 02:335, 02:429, 02:432, 02:450-02:452
*02:08
*02:09
*02:10, 02:244, 02:453
*02:11:01, 02:11:03-02:11:05, 02:297-02:298, 02:308
*02:11:02
*02:12
*02:13
*02:14
*02:15N
*02:16
*02:17:01-02:17:03, 02:303
*02:18
*02:19
*02:20:01-02:20:02
*02:21
*02:22:01-02:22:02
*02:24:01-02:24:02
*02:25
*02:26
*02:27
*02:28, 02:409 <sup>f</sup>
*02:29
*02:30
*02:31
*02:32N
Well No.
HLA-A allele <sup>1,2,3</sup>
Well No.
HLA-A allele <sup>1,2,3</sup>
*02:33
*02:34
*02:35:01
*02:35:02
*02:35:03

*02:36
*02:37
*02:38
*02:39
*02:40:01-02:40:02
*02:41
*02:42
*02:43N
*02:44
*02:45
*02:46
*02:47
*02:48
*02:49
*02:50
*02:51
*02:52
*02:53N
*02:54
*02:55
*02:56:01
*02:56:02
*02:57
*02:58
*02:59, 02:203
*02:60:01-02:60:02, 02:254
*02:61
*02:62
*02:63
*02:64
*02:65
*02:66
*02:67
*02:68
*02:69
*02:70
*02:71
*02:72
*02:73
*02:74:01-02:74:02
*02:75
*02:76:01-02:76:02
*02:77
Well No.
HLA-A allele <sup>1,2,3</sup>
Well No.
HLA-A allele <sup>1,2,3</sup>
*02:78
*02:79:01-02:79:02
*02:80
*02:81
*02:82N
*02:83N
*02:84
*02:85

<b>*02:86</b>
<b>*02:87</b>
*02:88N
*02:89
<b>*02:90, 02:331</b>
*02:91
*02:92
<b>*02:93</b>
*02:94N
*02:95
<b>*02:96</b>
<b>*02:97:01-02:97:02, 02:305N<sup>4</sup></b>
*02:99
<b>*02:101:01-02:101:02</b>
<b>*02:102</b>
*02:103
*02:104
<b>*02:105</b>
*02:106
<b>*02:107, 02:202, 02:251<sup>4</sup></b>
*02:108
*02:109
<b>*02:110</b>
*02:111, 02:350N <sup>4</sup>
*02:112
<b>*02:113N</b>
<b>*02:114</b>
*02:115
*02:116, 02:196 <sup>6</sup>
*02:117
<b>*02:118</b>
<b>*02:119, 02:263</b>
*02:120, 02:223N <sup>4</sup>
<b>*02:121, 02:425</b>
<b>*02:122</b>
<b>*02:123</b>
*02:124
<b>*02:125N, 02:217:02</b>
*02:126
*02:127
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>*02:128</b>
<b>*02:129</b>
*02:130
<b>*02:131</b>
<b>*02:132, 02:215, 02:237<sup>4</sup></b>
<b>*02:133</b>
*02:134, 02:314N <sup>4</sup>
*02:135
<b>*02:136</b>
*02:137



*02:138, 02:284N <sup>d</sup>
<b>*02:139, 02:235<sup>d</sup></b>
*02:140
<b>*02:141, 02:275<sup>d</sup></b>
*02:142
*02:143
*02:144
<b>*02:145, 02:221<sup>d</sup></b>
<b>*02:146</b>
<b>*02:147</b>
*02:148
<b>*02:149</b>
<b>*02:150, 02:197, 02:325<sup>d</sup></b>
<b>*02:151</b>
*02:152
<b>*02:153</b>
<b>*02:154</b>
*02:155
*02:156
*02:157
<b>*02:158</b>
<b>*02:159, 02:293Q, 02:364</b>
<b>*02:160</b>
<b>*02:161</b>
*02:162
<b>*02:163</b>
<b>*02:164</b>
<b>*02:165</b>
<b>*02:166</b>
<b>*02:167, 02:299</b>
*02:168
<b>*02:169</b>
*02:170
<b>*02:171:01-02:171:02</b>
<b>*02:172</b>
<b>*02:173</b>
<b>*02:174, 02:228</b>
*02:175
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*02:176
<b>*02:177</b>
*02:178
*02:179
<b>*02:180, 02:358<sup>d</sup></b>
*02:181
<b>*02:182</b>
<b>*02:183</b>
*02:184
*02:185
<b>*02:186</b>
*02:187
<b>*02:188</b>

*02:189
*02:190
*02:191
*02:192
*02:193, 02:213*
*02:194
*02:195
*02:198
*02:199
*02:200
*02:205
*02:206
*02:207
*02:211-01
*02:211-02
*02:214
*02:222N
*02:224, 02:371
*02:225N
*02:226N
*02:227N
*02:230
*02:233
*02:242
*02:245
*02:246
*02:247
*02:250N
*02:255
*02:259
*02:260
*02:267
*02:268
*02:269, 02:367
*02:271
Well No.
HLA-A allele <sup>1,2,3</sup>
Well No.
HLA-A allele <sup>1,2,3</sup>
*02:274
*02:279
*02:280
*02:282
*02:283
*02:289
*02:291
*02:295
*02:300
*02:301N
*02:304
*02:309
*02:310
*02:315
*02:319
*02:321N

*02:322
<b>*02:323</b>
*02:324
<b>*02:330</b>
*02:334
<b>*02:338</b>
*02:339
<b>*02:342</b>
<b>*02:344</b>
<b>*02:351</b>
*02:355
*02:356N
*02:366N
*02:369
*02:373N
*02:376
<b>*02:377</b>
*02:382
*02:387
<b>*02:393</b>
*02:395N
*02:398
*02:400
*02:402
*02:403
*02:406
<b>*02:407</b>
*02:408
*02:412
*02:413
*02:415
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*02:417
<b>*02:419</b>
*02:420
*02:421
<b>*02:426</b>
*02:433
*02:436
<b>*02:437</b>
*02:438
*02:439N
*02:440Q
*02:447
*02:449
*02:454

*01:13, 01:17, 03:01:01-01-03:01:02, 03:01:04-03:07, 03:10-03:11N, 03:13-03:16, 03:17:02-03:22:02, 03:23:02-03:31, 03:33-03:35, 03:37-03:40, 03:42-03:43, 03:45-03:50, 03:52-03:58, 03:60-03:71, 03:73-03:74, 03:76-03:87, 03:90-03:94, 03:96-03:106, 03:109-03:122, 03:124-03:141, 03:143-03:151, 03:153-03:169, 03:174-03:175, 34:08, 74:13
*01:20, 01:66, 01:130
*01:44
*01:69:02
*01:83:02, 11:121, 26:07:02, 33:01:02, 33:08-33:09, 74:04
*01:104
*01:129, 29:15, 31:78
*01:134, 03:108, 03:172
*03:01:03
*03:09
*03:17:01
*03:23:01
*03:44
*03:51
*03:75
*03:89
*03:95
*03:123:01-03:123:02
*03:170
*03:171
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*11:06, 11:18
*11:16, 11:35, 11:57, 11:73, 11:158, 29:51
*11:119:01-11:119:02
*23:01:01-23:01:11, 23:05-23:07N, 23:09, 23:11N, 23:13, 23:14:02-23:23, 23:25-23:30, 23:32-23:33, 23:35-23:40, 23:42, 23:44, 23:46, 23:48-23:51, 23:53-23:54, 23:56, 23:58-23:59, 24:24
*23:02
*23:03:01-23:03:02, 29:07, 31:29
*23:04, 23:43
*23:08N
*23:10

*23:12
*23:14:01
*23:24, 23:57, 24:13:01
*23:31, 23:45, 23:55
*23:34, 24:02:16, 24:13:02
*23:41
*23:47
*23:52
*23:60
*24:02:01:01-24:02:15, 24:02:17-24:02:41, 24:02:43- 24:02:69, 24:02:71, 24:04- 24:05, 24:07-24:09N, 24:11N, 24:17, 24:19-24:20, 24:25- 24:27, 24:29, 24:31-24:32, 24:34-24:40N, 24:43-24:45N, 24:47-24:50, 24:52-24:53, 24:57-24:58, 24:60N-24:64, 24:66-24:70, 24:74:01- 24:74:02, 24:76-24:86N, 24:88, 24:90N, 24:95-24:98, 24:100-24:106, 24:108-24:109, 24:111-24:124, 24:126-24:130, 24:132N, 24:134-24:137, 24:139-24:156, 24:159, 24:161 24:166, 24:168-24:176, 24:178 24:187, 24:189, 24:191- 24:203, 24:205-24:206, 24:209, 24:212, 24:214- 24:217, 24:220-24:226, 24:230 24:233, 24:236-24:252N, 33:19
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*24:02:70, 24:03:01-24:03:02, 24:33, 24:75
*24:06, 24:54, 24:87, 24:107, 24:160, 24:213, 24:227
*24:10:01-24:10:02
*24:14, 24:93
*24:15, 24:41, 24:51, 24:92, 24:235
*24:18
*24:23
*24:28, 24:30, 24:42, 24:89
*24:46
*24:55
*24:56, 24:177
*24:59, 24:190, 24:229
*24:71
*24:72

*24:73, 24:157
*24:91
*24:94
*24:99
*24:131, 24:218
*24:133
*24:138
*24:188
*24:207
*24:210
*24:219
*24:228
*24:234
*25:11, 26:36
*26:03:01-26:03:02, 26:06, 26:21, 26:30, 26:78, 26:92
*26:10
*26:22, 66:09
*29:19, 32:08, 33:24
*29:22, 33:22
*29:48
*30:01:01-30:01:07, 30:08, 30:11:01-30:11:02, 30:14L- 30:15, 30:17-30:20, 30:23- 30:24, 30:26, 30:30-30:31, 30:35-30:36, 30:38-30:43, 30:48-30:49, 30:53-30:54, 30:56, 30:58-30:60, 30:62- 30:63, 30:65, 30:72-30:73N, 30:75, 30:78N, C*03:82
*30:13, 30:44, 30:46
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*30:16
*30:22
*30:37
*30:52
*30:55
*30:61
*30:74
*31:48, 33:15
*31:67-31:68, 33:32:01
*32:28
*32:46
*32:48N

*68:01:01:01-68:01:10, 68:01:12-68:01:13, 68:01:15- 68:01:20, 68:03:01-68:04, 68:07, 68:10-68:11N, 68:13:01 68:13:02, 68:16-68:17, 68:19, 68:21:01-68:22, 68:24, 68:29, 68:32, 68:35-68:39, 68:41, 68:43, 68:45-68:47, 68:52, 68:55-68:59N, 68:66, 68:68- 68:70, 68:72-68:73, 68:75:01- 68:75:02, 68:79, 68:84-68:85, 68:87-68:91, 68:93-68:95, 68:98-68:104, 68:106
*68:01:11
*68:01:14
*68:02:01:01-68:02:05, 68:18N, 68:27:01-68:28, 68:31, 68:34, 68:40, 68:53, 68:60, 68:62, 68:64, 68:67, 68:74, 68:77-68:78, 68:80, 68:82, 68:86, 68:92, 68:97
*68:05, 68:20
*68:06
*68:08:01-68:08:02
*68:09
*68:12
*68:14
*68:15
*68:23
*68:25
*68:26
*68:30
*68:33, 68:96
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
<b>Well No.</b>
<b>HLA-A allele<sup>1,2,3</sup></b>
*68:42
*68:44
*68:48
*68:49N
*68:50
*68:51
*68:54
*68:61
*68:63
*68:65
*68:71
*68:76
*68:81
*68:83
*68:105
*69:01-69:02
*74:06

*80:01:01:01
B*15:67, B*35:110
B*15:173, B*18:63
B*44:59:02, B*44:136
B*51:136
C*07:204
C*12:37
C*12:94, C*14:48
<b>Well No.</b>
HLA-A allele <sup>1,2,3</sup>