

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32						
HLA-A allele <sup>1,2,3</sup>																																						
*23:01:01-23:01:13	1	2																18			21																	
*23:02		2											13					w			?									?								
*23:03:01-23:03:02	1		3																		?										?							
*23:04	1												14					18			?										?							
*23:05	1	2							9									18											28									
*23:06	1	2			5													18			21																	
*23:07N	1	2			6													18			21																	
*23:08N, 23:22 <sup>4</sup>	1	2					7											18			?										?							
*23:09, 23:26 <sup>4</sup>	1	2						8										18			21																	
*23:10, 23:23 <sup>4</sup>	1	2								10								18			?											?						
*23:11N	1	2									11							18			?												?					
*23:12	1	2										12						18			?												?					
*23:13, 23:33 <sup>4</sup>	1	2															17	18		?													?					
*23:14:01	1	2													15	16		w		?													?					
*23:14:02	1	2														16		18		?														?				
*23:15, 23:27 <sup>4</sup>	1	2																18	19	?														?				
*23:16, 23:29 <sup>4</sup>	1	2																18		20	?													?				
*23:17	1	2																18												28								
*23:18, 23:28 <sup>4</sup>	1	2		4														18		?														?				
*23:19Q, 23:31, 23:53 <sup>4</sup>	1	2																18		?	22													?				
*23:20	1	2																18		?		23												?				
*23:21	1	2																18		?			24											?				
*23:24	1	2											13					18		?														?				
*23:25	1	2	3															18		?															?			
*23:30	1	2											13					18		?															?			
*23:32	1	2									11							18		20	?														?			
*23:34, 23:57	1	2																18		?															27	?		
*23:35	1	2																18		?														27	?			
*23:36	1	2																18		?								26					?					
*23:37:01-23:37:02	1	2																18		?									?			29						
*23:38N	1	2																18		21															30			
*23:39	1	2					7	8										18		?															?			
*23:40	1	2																18		?			24										?	29				
*23:41-23:42 <sup>4</sup>	1	2																18		?				25									?					
*23:43	1									10								18		?															?			
*23:44	1	2																18		?																31		
*23:45	1	2																18		?	22	23													?			
*23:46	1	2																18	19	20	?														?			
*23:47	1	2													15			18		?														?	29			
*23:48	1	2						8										18		?														?	29			
*23:49	1	2																18	19	?			24											?				
*23:50-23:52, 23:54-23:55, 23:59-23:60	1	2																18		?														?				
*23:56	1	2																18		?														?		30		
*23:58	1	2		4														18		?			23											?				
*01:02, 01:20								8																														
*01:04N, 03:21N, 11:21N					6																																	
*02:17:01-02:17:03	1													14		w																27						
*02:19, 02:36-02:37, 02:54, 02:255, 02:417, 68:26, 68:65		2																																				
*02:40:01, 02:51, 02:130																		18																				
*02:40:02, 68:51																		w																				
*02:41, 02:80, 02:117, 02:289, 02:304, 02:454, 26:10																							23															
*02:82N							7																															
*02:108, 02:110, 02:268, 02:300, 02:303	1												14																							27		
*02:221, 29:15, 31:78																											25											
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32						

Negative Control



*24:46		2							10				16																								?						
*24:61		2											16																							28		30					
*24:62		2											16									22	23																?				
*24:66		2											16		19	20																							?				
*24:77		2	4										16																											?			
*24:94	1												14																								27	?					
*24:125, 24:204													14	16																										?			
*24:128		2											16			20																								28			
*24:129		2					8						16																											28			
*24:138													13	14																											?		
*24:167													13	14	16																											?	
*24:188, 24:228	1	2																																							27	?	
*24:207	1												14	16																											27	?	
*24:208			3											16								22																			?		
*24:234		2											15	16																												28	29
*24:255		2											16		19																											28	
*24:267		2	4										16																													28	
*25:11, 32:08										12																																	
*26:16														16																													
*29:07, 29:49, 31:29	1	3											14									22																					
*30:06																19							24																				
*30:22					7	8																																					
*30:75															17	19	20																										
*31:36					5																																						
*31:67-31:68, 32:28, 33:32:01																18							23																				
*32:05																																											26
*68:45														16	17																												
B*18:27	1	3							10				14	15																													
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32											

Negative Control