



HLA-A23 2018 0G9

*23:36	1	2															18		?					26	?														
*23:37:01-23:37:02	1	2															18		?						?	29													
*23:38N	1	2															18		21								30												
*23:39	1	2					7	8									18		?							?													
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>	<b>25</b>	<b>26</b>	<b>27</b>	<b>28</b>	<b>29</b>	<b>30</b>	<b>31</b>	<b>32</b>	<b>33</b>	<b>34</b>	<b>35</b>	<b>36</b>			
<b>Well No.</b>	<b>1</b>	<b>2</b>	<b>3</b>	<b>4</b>	<b>5</b>	<b>6</b>	<b>7</b>	<b>8</b>	<b>9</b>	<b>10</b>	<b>11</b>	<b>12</b>	<b>13</b>	<b>14</b>	<b>15</b>	<b>16</b>	<b>17</b>	<b>18</b>	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>	<b>25</b>	<b>26</b>	<b>27</b>	<b>28</b>	<b>29</b>	<b>30</b>	<b>31</b>	<b>32</b>	<b>33</b>	<b>34</b>	<b>35</b>	<b>36</b>			
*23:40	1	2															18		?				24			?	29												
*23:41-23:42 <sup>4</sup>	1	2															18		?				25			?													
*23:43	1								10								18		?								?												
*23:44	1	2															18		?								?			31									
*23:45	1	2															18		?	22	23						?												
*23:46	1	2															18	19	20	?							?												
*23:47	1	2												15			18		?								?	29											
*23:48	1	2						8									18		?								?	29											
*23:49	1	2															18		?			24				?										35			
*23:50-23:52, 23:54-23:55, 23:60-23:64, 23:67, 23:71-23:75, 23:81-23:82	1	2															18		?								?												
*23:53	1	2											14				18		?	22						27	?												
*23:56	1	2															18		?								?			30									
*23:58	1	2															18		?		23					?										34			
*23:59	1	2															18		?								?								32				
*23:65	1	2															18		21																		33		
*23:66		2											13						?								?												
*23:68	1	2												14			18		21																				
*23:70	1		3											14					?	22						27	?												
*23:79	1																18		21																				
*23:83	1		3											14					?								?												
*01:02, 01:20, 01:188, 01:190								8																															
*01:04N, 03:21N, 11:21N						6																																	
*01:88:02, 01:236, 02:19, 02:36-02:37, 02:54, 02:255, 02:417, 11:209, 68:26, 68:65, 68:115, 68:131		2																																					
*01:203, 02:671, 11:166, 31:85, 32:72, 80:01:01:01-80:03																																						32	

Negative Control









HLA-A23 2018 0G9

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36		
*30:22							7	8																														
*30:56																	17																					
*30:75																	17		19	20																		
*30:104							7	8									17																					
*31:29	1		3											14								22																
*31:36					5																																	
*31:67-31:68, 32:28, 32:66, 33:32:01																		18				23																
*32:05, 32:79																										26												
*68:45, 68:117																16	17																					
*68:51																		w																				
B*18:27	1		3							10				14	15																							
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36		

Negative Control