

INTERPRETATION TABLE

HLA-A*23 Add-on SSP typing

	Well		
	1	2	3
Length of spec.	210	90	110
PCR product			
Length of int.	800	1070	1070
pos. control ¹			
5'-primer ²	368	920	453
	5' -gTT 3'	5' -CCA 3'	5' -AAA 3'
3'-primer ³	539	971	524
	5' -TCA 3'	5' -CAg 3'	5' -CAC 3'
Well No.	1	2	3
HLA-A allele			
*23:01:01-23:01:07, 23:06-23:07N, 23:09, 23:26, 23:38N	1	2	3
*23:02		?	w
*23:03:01-23:03:02, 23:24, 23:34	1	?	
*23:04, 23:08N, 23:10-23:13, 23:14:02-23:16, 23:18-23:23, 23:25, 23:27-23:33, 23:35-23:37:02, 23:39-23:50	1	?	3
*23:05, 23:17, 24:24	1		3
*23:14	1	?	w
*02:17:01-02:17:02, 02:108, 02:110, 02:268, 02:300, 02:303, 24:13:01-24:13:02, 24:18, 24:94, 24:188, 29:07, 31:29, B*18:27	1		
*02:40, 02:51, 02:130, 32:28, 33:32			3
*68:51			w
HLA-A allele			
Well No.	1	2	3