

<b>INTERPRETATION TABLE</b>			
<b>HLA-A*23 Add-on SSP typing</b>			
	<b>Well</b>		
	<b>1</b>	<b>2</b>	<b>3</b>
Length of spec.	210	90	110
PCR product		235	
Length of int.	800	1070	1070
pos. control <sup>1</sup>			
5'-primer <sup>2</sup>	368	28	453
	5' -gTT 3'	5' -TCg 3'	5' -AAA 3'
		920	
		5' -CCA 3'	
3'-primer <sup>3</sup>	539	92	524
	5' -TCA 3'	5' -AAC 3'	5' -CAC 3'
		971	
		5' -CAg 3'	
Well No.	1	2	3
HLA-A allele			
*23:01:01-23:01:11, 23:05-23:07N, 23:09, 23:26, 23:38N	1	2	3
*23:02		?	W
*23:03:01-23:03:02, 23:24, 23:34, 23:57	1	?	
*23:04, 23:08N, 23:10-23:13, 23:14:02-23:16, 23:18-23:23, 23:25, 23:27-23:33, 23:35-23:37:02, 23:39-23:56, 23:58	1	?	3
*23:14:01	1	?	W
*23:17, 24:24	1		3
*24:25		2	
*02:17:01-02:17:03, 02:108, 02:110, 02:268, 02:300, 02:303, 24:13:01-24:13:02, 24:18, 24:94, 24:188, 24:207, 24:228, 29:07, 29:49, 31:29, B*18:27	1		
*02:40:01, 02:51, 02:130, 31:67-31:68, 32:28, 33:32:01			3
*02:40:02, 68:51			W
HLA-A allele			
Well No.	1	2	3