

INTERPRETATION TABLE		
HLA-A*2409N SSP typing		
	Well	
	1	2
Length of spec.	105	175
PCR product		205
Length of int.	800	1070
pos. control¹		
5'-primer(s)²	678	98
	5'-AgA ^{3'}	5'-CTC ^{3'}
		368
		5'-gTT ^{3'}
3'-primer(s)³	742	259
	5'-CTA ^{3'}	5'-gTT ^{3'}
		502
		5'-CTT ^{3'}
		539
		5'-TCT ^{3'}
Well No.	1	2
HLA-A allele		
*2409N	1	2
*021701, 021702		w
*2314, 24020101-240217, 240220-2408, 2410, 2411N, 241301, 241302, 2417-2450, 2454-2456, 2458-2463, 2466- 2491, 2493, 2495-2499, 2616, 3319, 6845		2
HLA-A allele		
Well No.	1	2