

INTERPRETATION TABLE

HLA-A*24:09N SSP typing

	Well	
	1	2
Length of spec.	105	175
PCR product		205
Length of int.	800	1070
pos. control ¹		
5'-primer(s) ²	678	98
	5'-AgA 3'	5'-CTC 3'
		368
		5'-gTT 3'
3'-primer(s) ³	742	259
	5'-CTA 3'	5'-gTT 3'
		502
		5'-CTT 3'
		539
		5'-TCT 3'
Well No.	1	2
HLA-A allele		
*24:09N	1	2
*02:17:01-02:17:02		w
*23:14, 24:02:01:01-24:08, 24:10-24:11N, 24:13:01-24:13:02, 24:17-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139-24:182, 26:16, 33:19, 68:45		2
HLA-A allele		
Well No.	1	2