

## INTERPRETATION TABLE

### HLA-A\*24:09N SSP typing

	Well	
	1	2
Length of spec.	105	175
PCR product		205
Length of int.	<b>800</b>	1070
pos. control <sup>1</sup>		
5'-primer(s) <sup>2</sup>	<b>678</b>	<b>98</b>
	5' -AgA 3'	5' -CTC 3'
		<b>368</b>
		5' -gTT 3'
3'-primer(s) <sup>3</sup>	<b>742</b>	<b>259</b>
	5' -CTA 3'	5' -gTT 3'
		<b>502</b>
		5' -CTT 3'
		<b>539</b>
		5' -TCT 3'
Well No.	<b>1</b>	<b>2</b>
HLA-A allele		
<b>*24:09N</b>	<b>1</b>	<b>2</b>
<b>*02:17:01-02:17:02</b>		<b>w</b>
<b>*23:14:01-23:14:02, 24:02:01:01-24:08, 24:10-24:11N, 24:13:01-24:13:02, 24:17-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139- 24:187, 24:189-24:209, 26:16, 33:19, 68:45</b>		<b>2</b>
HLA-A allele		
Well No.	<b>1</b>	<b>2</b>