

## INTERPRETATION TABLE

### HLA-A\*24:09N SSP typing

	Well	
	1	2
Length of spec.	105	175
PCR product		205
Length of int.	<b>800</b>	1070
pos. control <sup>1</sup>		
5'-primer(s) <sup>2</sup>	678	98
	5' -AgA 3'	5' -CTC 3'
		368
		5' -gTT 3'
3'-primer(s) <sup>3</sup>	742	259
	5' -CTA 3'	5' -gTT 3'
		502
		5' -CTT 3'
		539
		5' -TCT 3'
Well No.	1	2
HLA-A allele		
*24:09N	1	2
*02:17:01-02:17:02		w
*23:14:01-23:14:02, 24:02:01:01-24:08, 24:10:01-24:11N, 24:13:01-24:13:02, 24:17-24:50, 24:54-24:56, 24:58-24:63, 24:66-24:91, 24:93, 24:95-24:113, 24:115-24:137, 24:139- 24:187, 24:189-24:210, 24:212-24:214, 24:216-24:221, 26:16, 33:19, 68:45		2
HLA-A allele		
Well No.	1	2