

HLA-A25 2019 5H3

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
HLA-A allele^{1,2,3}																
*25:01:01:01 -25:01:01:05, 25:01:03 -25:01:14, 25:14, 25:18, 25:20-25:23, 25:25-25:26, 25:28-25:29, 25:31-25:32, 25:34-25:35, 25:37-25:41, 25:45, 25:47-25:48, 25:51, 25:53	1	2	3											14		
*25:01:02, 25:33	1	2	3													
*25:02	1	2		4	5									14		
*25:03	1	2	3			6								14		
*25:04, 25:44	1	2	3				7							14		
*25:05	1	2	3					8						14		
*25:06		2	3				7							14		
*25:07	1	2	3							10				14		
*25:08	1	2	3						9					14		
*25:09	1	2	3								11			14		
*25:10	1	2	3									12		14		
*25:11	1	2	3										13	14		
*25:12N	1	2	3				7	8						14		
*25:13	1	2													15	
*25:15, 25:17	1													14		
*25:16	1	2	3					8	9					14		
*25:19:01 , 25:30	1	2	3			6	7							14		
*25:19:02	1		3			6	7							14		
*25:24	1	2	3	4										14		
*25:27:01-25:27:02	1	2	3					8		10				14		
*25:36, 25:43, 25:50, 25:52, 32:62 ⁴		2	3											14		
*25:42N	1	2	3				7		9					14		
*25:46	1		3											14		
*25:49N	1	2	3						9	10				14		
*01:01:01:01-01:01:27, 01:01:29-01:01:56, 01:01:58-01:01:61, 01:01:63-01:04:01:02N, 01:06, 01:08-01:12, 01:14-01:16N, 01:18N-01:33, 01:35-01:50, 01:52:01N-01:70, 01:72-01:99, 01:101-01:104, 01:106-01:143, 01:145-01:175, 01:177-01:189, 01:191-01:193, 01:195-01:203, 01:205-01:249, 01:251-01:281, 02:346, 02:427, 02:749, 03:41, 11:17, 11:19, 11:40, 11:98, 11:121, 26:04, 26:29, 26:49, 26:75, 26:78, 33:13, 36:01-36:08, 74:10, 74:30, 80:01:01:01-80:04			3													
*01:13, 01:17, 01:176, 03:63, 03:88, 11:01:01:01-11:01:67, 11:01:69-11:01:72, 11:01:74-11:01:75, 11:01:77-11:11, 11:13-11:16, 11:20-11:27, 11:29-11:39, 11:41-11:52Q, 11:54:01-11:95, 11:97, 11:99N-11:105, 11:107-11:115N, 11:117-11:120, 11:122-11:139, 11:141-11:158, 11:160-11:177, 11:179-11:249, 11:251N-11:290, 11:292-11:300, 11:302-11:306, 29:66, 30:125, 30:140, C*07:335				4												
*01:51, 02:644			3			6										

Negative Control

*32:15		2			6													
*32:26:01-32:26:02	1	2																
*34:01:01-34:01:02, 34:05, 34:12, 34:16				4	5		7											
*34:02:01-34:02:04, 34:04, 34:08, 34:13, 34:15				4	5	6	7					12						
*34:03				4	5	6	7					11	12					
*34:06				4	5	6	7					11						
*34:07					5	6	7						12					
*34:09					5	6	7						12					15
*34:10N				4	5								12					
*34:11, 66:10, 69:02				4	5													
*34:14				4			7											
*34:17				4	5		7					11						
*66:02, 66:16, 66:21, 66:25-66:26Q					5													15
*66:03, 66:28N																		15
*66:06	1			4	5	6	7											
*66:12	1				5													15
*66:14	1			4	5			8		10								
*66:18	1			4	5					9								
*66:22	1			4	5							11						
*68:71					5	6												
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	Negative Control	