

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48						
HLA-A allele <sup>1,2,3</sup>																																																						
*26:01:01:01-26:01:01:02, 26:01:01:04-26:01:20, 26:01:22-26:01:39, 26:01:41, 26:77, 26:79-26:80, 26:83-26:85, 26:87-26:90, 26:93-26:97, 26:100-26:102, 26:105, 26:108-26:110, 26:113, 26:117, 26:119-26:123, 26:125-26:126, 26:128-26:129, 26:131	1	2																																																				
*26:01:01:03N	1	2																																															47					
*26:01:21	1	2																								26																												
*26:02:01-26:02:02, 26:118, 26:132			2	3																																																		
*26:03:01, 26:92, 26:111	1			4																									28																									
*26:04		2			5																																																	
*26:05	1					6																																																
*26:06				4								11																	28																									
*26:07:01-26:07:02	1										10																																											
*26:08, 26:67, 26:81, 26:86, 26:124	1	2					7																																															
*26:09		2																			17																																	
*26:10	1	2												14																																								
*26:11N	1	2						8																																														
*26:12	1	2							9																																													
*26:13	1	2													15	16																																						
*26:14	1	2						8				12																																										
*26:15	1	2																																																				
*26:16	1	2										12																																										
*26:17, 26:106 <sup>4</sup>	1	2																																																				
*26:18	1	2						8	9		12			13																																						w		
*26:19	1										11				15	16																																						
*26:20	1	2								10																																												
*26:21	1			4															18										28																	32								
*26:22	1	2																		19																																		
*26:23, 26:27 <sup>4</sup>	1	2																				20																																
*26:24, 26:41 <sup>4</sup>	1	2																						21																														
*26:25N, 26:38 <sup>4</sup>	1	2																						22																														
*26:26	1	2																							23																													
*26:28	1	2					7	8				12																																								46		
*26:29	w	2								9																																												
*26:30	1			4																																																		
*26:31	1																	18										25																										
*26:32, 26:70	1	2																																																				
*26:33	1	2													15	16																																						
*26:34		2																																																				
*26:35	1	2																																																				
*26:36, 26:116 <sup>4</sup>	1	2																																																				
*26:37	1	2																																																				
*26:39	1	2																																																				
*26:40		2																	19																																			
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48						
Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48						
*26:42	1	2																																																				
*26:43:01, 26:61 <sup>4</sup>	1	2																																																				
*26:43:02	1	2																																																				
*26:45	1	2										11	13																																									

Negative Control







