

HLA-A29 2017 9F9

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28
HLA-A allele^{1,2,3}																												
*29:01:01:01 , 29:01:02-29:01:06, 29:01:08-29:01:09, 29:28, 29:57-29:58, 29:60-29:62, 29:67, 29:76, 29:79, 29:83, 29:99, 29:101	1													14	15													
*29:01:01:02N	1	2												14	15													
*29:01:07														14	15													
*29:02:01:01-29:02:03 , 29:02:05-29:02:21, 29:30-29:31, 29:38, 29:41-29:42, 29:45, 29:47, 29:50, 29:52, 29:59, 29:63, 29:65-29:66, 29:68, 29:70, 29:72, 29:74-29:75, 29:80, 29:84-29:86, 29:88-29:89, 29:91, 29:93-29:97, 29:100, 29:102			3											14	15													
*29:02:04												12		14	15													
*29:03			3	4										14	15													
*29:04			3		5									14	15													
*29:05, 29:87			3			6								14	15													
*29:06, 29:34⁴			3				7							14	15													
*29:07, 29:46⁴			3					8						14	15													
*29:08N			3						9					14	15													
*29:09			3							10				14	15													
*29:10:01-29:10:02			3								11			14	15													
*29:11, 29:92⁴			3										13	14	15													
*29:12	1												13	14	15													
*29:13			3											14							22							
*29:14			3								11				15													
*29:15	1													14	15	16												
*29:16	1								9					14	15													
*29:17	1									10				14	15													
*29:18												12			15													
*29:19			3			6	7							14														
*29:20, 29:82⁴	1						7							14	15													
*29:21			3											14	15	16												
*29:22			3											14	15						21							
*29:23			3		5						11			14	15													
*29:24	1													14	15		17											
*29:25														14	15				19									
*29:26			3											14	15				19									

tive Control

