

HLA-A29 2020 5L0

Well No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
HLA-A allele^{1,2,3}																															
*29:01:01:01 , 29:01:01:03-29:01:06, 29:01:08-29:01:10, 29:01:12, 29:28, 29:57-29:58, 29:60-29:62, 29:67, 29:76, 29:79, 29:83, 29:99, 29:101:01-29:101:02, 29:107, 29:110, 29:113, 29:115, 29:117-29:118, 29:124-29:125	1													14	15																
*29:01:01:02N	1	2												14	15																
*29:01:07, 29:01:13, 29:109, 29:136														14	15																
*29:01:11	1											12		14	15																
*29:02:01:01-29:02:03 , 29:02:05-29:02:29, 29:30-29:31, 29:38, 29:41-29:42, 29:45, 29:47, 29:50, 29:52, 29:59, 29:63, 29:65-29:66, 29:68, 29:70, 29:72, 29:74-29:75, 29:80, 29:84-29:86, 29:89, 29:91, 29:93-29:97, 29:100, 29:102, 29:106, 29:108, 29:111, 29:114, 29:116, 29:119-29:122, 29:126Q-29:132, 29:134-29:135, 29:137-29:138, 29:140-29:143			3											14	15																
*29:02:04												12		14	15																
*29:03			3	4										14	15																
*29:04			3		5									14	15																
*29:05, 29:87 , 29:104			3			6								14	15																
*29:06			3				7							14	15																
*29:07, 29:46 ⁴			3					8						14	15																
*29:08N, 29:139 ⁴			3						9					14	15																
*29:09			3							10				14	15																
*29:10:01-29:10:02			3								11			14	15																
*29:11 , 29:92 ⁴			3										13	14	15																
*29:12	1												13	14	15																
*29:13			3											14							22										
*29:14			3								11			15																	
*29:15	1													14	15	16															
*29:16	1								9					14	15																
*29:17, 29:144 ⁴	1									10				14	15																
*29:18												12		15																	
*29:19			3											14															29		
*29:20, 29:82⁴	1													14	15													29			
*29:21			3											14	15	16															
*29:22			3											14	15						21										
*29:23			3		5						11			14	15																

Control

